CT-U.S.-EASIA

KAP0421

Read entire protocol before use.

CT-U.S.-EASIA

I. INTENDED USE

Immunoenzymetric assay for the *in vitro* quantitative measurement of human Calcitonin (CT) in serum

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource CT-U.S.-EASIA Kit
- B. **Catalogue number :** KAP0421 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A Biological activities

Calcitonin(CT) is a 32 amino acid peptide hormone secreted by the para-follicular C-cells of the thyroid gland under serum calcium control. After acute administration, this peptide acts as a potent hypocalcemic and hypophosphatemic hormone by increasing renal calcium clearance and reducing bone resorption. However, its precise physiological role in bone metabolism is not yet fully understood.

Various forms of CT may be detected in blood samples, including a CT monomer, an oxidized monomer, a dimer, higher molecular weight forms, and possibly precursor of CT. The concentrations of these peptides vary with clinical status, renal function and tissular origin of CT (normal or ectopic production).

Medullary thyroid carcinoma (MTC) is a malignant tumor, developed from the C-cells, secreting calcitonin in large excess. This disease occurs either as a sporadic (80%) or a familial (20%) form, which is transmitted as an autosomal dominant gene or as a component of multiple endocrine neoplasia (IIB).

Moderate hypercalcitoninemia is also observed in pregnancy, pernicious anaemia, renal failure and during the neonatal period. Preferably, monomer form of CT is detected in this assay.

B Clinical Application

The measurement of CT is used for :

- Diagnosis of medullary thyroid carcinoma (MTC),
- Follow up of malignant tumors, to check the success of surgery and to monitor for recurrence,
- Diagnosis of the preclinical cases of the familial forms of MTC (MEN II or Sipple syndrome) by the use of stimulation tests (calcium or pentagastrin),
- Study of the pathophysiology of the calcium-phosphate and bone metabolism.

IV. PRINCIPLES OF THE METHOD

The DIASource CT-U.S.-EASIA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplate. Calibrators and samples react with the capture monoclonal antibody (MAb 1) coated on microtiter well and with a monoclonal antibody (MAb 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAb 1 – human CT – MAb 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB ready for use) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the Calcitonin concentration.

A calibration curve is plotted and Calcitonin concentration in samples is determined by interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Color Code	Reconstitution
 Microtiterplate with 96 anti CT (monoclonal antibodies) coated breakable wells	96 wells	blue	Ready for use
AbHRPCONC Conjugate: HRP labelled anti-CT (monoclonal antibodies) in Stabilizing Buffer	1 vial 0.125 ml	yellow	Dilute 50 x with conjugate buffer
CONJBUF Conjugate buffer: TRIS-Maleate buffer with bovine serum albumin, EDTA and thymol	1 vial 6 ml	red	Ready for use
CALN Calibrator N = 0 to 5 (see exact values on vial labels) in CT-Free human serum	6 vials lyophil.	yellow	Add 0.5 ml distilled water
SERUM CT free human serum (to be used for samples dilution) with thymol	1 vial lyophil.	black	Add buffer (see reconstitution volume on the label)
BUF Buffer (serum free): borate buffer	1 vial 8 ml	black	Ready for use
WASHSOLNCONC Wash Solution (Tris-HCl)	1 vial 10 ml	brown	Dilute 200 x with distilled water (use a magnetic stirrer).
CONTROLN Controls - N = 1 or 2 in human serum with gentamycin	2 vials lyophil.	silver	Add 0.5 ml distilled water
CHROMTMB Chromogenic TMB solution (Tetramethylbenzidine)	1 vial 12 ml	brown	Ready for use
STOPSOLN Stop Solution: HCl: 1N	1 vial 12 ml	white	Ready for use

Note: 1. CT free human serum is to be used for samples dilution.
2. 1 pg of our reference preparation is equivalent to 0.19 µIU NIBSC 89/620.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 50 µl, 100 µl, 500 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Washer for microtiterplate
- Microtiterplate reader capable of reading at 450 nm and 650 nm (bichromatic reading).

VII. REAGENT PREPARATION

- Calibrators** : Reconstitute the calibrators with 0.5 ml distilled water.
- Controls** : Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.
- CT Free Serum**: Reconstitute the CT Free Serum with the amount of Buffer as indicated on the vial label. Allow it to remain undisturbed until completely dissolved, and then mix well by gentle inversion.
- Working anti-CT-HRP conjugate**: Prepare an adequate volume of conjugate solution by adding for example :40 µl of the 50 x concentrated anti-CT-HRP conjugate to 2 ml of conjugate buffer. Use a vortex to homogenize. Extemporaneous preparation is recommended.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused wells must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, calibrators, controls and CT free serum should be frozen immediately after use and kept at –20°C for 3 months. Only one freeze thawing cycle is allowed, discard the calibrators, controls and CT free serum after the second use.
- The concentrated Wash Solution is stable at room temperature until expiration date.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, the concentrated conjugate (50x) is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- The Working anti-CT-HRP conjugate is stable for 1 week at 4°C.
- The chromogenic TMB solution and the Stop Solution are stable at 2°C to 8°C until the expiry date.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum must be kept at 2 - 8°C.
- If the test is not run within 24 hours, storage in aliquots at -20°C is recommended. Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.
- Do not use haemolysed samples.
- Do not use lipemic samples.

X. PROCEDURE

- Handling notes**
Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.
High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
Prepare a calibration curve for each run, do not use data from previous runs.
The chromogenic solution should be colourless. If a dark blue colour develops within a few minutes after preparation, this indicates that the preparation is unusable and must be discarded.

Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.

During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

B. Procedure

1. Select the required number of wells for the run. The unused wells should be resealed in the bag with a desiccant and stored at 2-8°C.
2. Secure the wells into the holding frame.
3. Pipette 100 µl of each Calibrator, Control and Sample into the appropriate wells.
4. Pipette 50 µl of Working anti-CT-HRP conjugate into all the wells.
5. Incubate for 18 ± 1 hour at 2-8°C.
6. Aspirate the liquid from each well.
7. Wash the plate 3 times by:
 - Dispensing 0.4 ml of Wash Solution into each well
 - Aspirating the content of each well
8. Pipette 100 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
9. Incubate the microtiterplate for 30 minutes at room temperature avoiding direct sunlight.
10. Pipette 100 µl of Stop Solution into each well.
11. Read the absorbances at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section XI.

XI. CALCULATION OF RESULTS

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of Calcitonin (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
4. Read the concentration for each control and sample by interpolation on the calibration curve.
5. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

CT-U.S.-EASIA		Bichromatic model (OD)
Calibrator	0 pg/ml	0.009
	10 pg/ml	0.029
	50 pg/ml	0.127
	100 pg/ml	0.447
	200 pg/ml	0.919
	400 pg/ml	1.87

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 0.7 pg/ml.

B. Specificity

Some potentially interfering hormones have been tested in this assay. At concentrations up to 100 ng/ml, none of the following hormones showed significant interference :

- CGRP
- Salmon-calcitonin
- PDN 21
- Procalcitonin N terminal.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> ± SD (pg/ml)	CV (%)	Serum	N	<X> ± SD (pg/ml)	CV (%)
A	19	43.0 ± 0.75	1.7	A	8	44.6 ± 2.1	4.9
B	19	133.7 ± 5.2	3.9	B	8	136.3 ± 8.1	6

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Added CT (pg/ml)	Recovered CT (pg/ml)	Recovery (%)
327.7	340.6	104
160.7	159.3	99
80.5	80.4	99
48.3	50.8	105

DILUTION TEST

Sample	Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)
Serum 1	1/1	-	300.6
	1/2	150.3	157.9
	1/4	75.1	75.5
	1/8	37.6	45.7
	1/16	18.8	25.2
	1/32	9.4	12.1
	1/64	4.7	5

Samples were diluted with CT Free Serum.

E. Hook effect

A sample spiked with CT up to 480000 pg/ml gives higher OD's than the last calibrator point.

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls that contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Normal values

84 samples from normal subjects obtained values below 11 pg/ml.

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl, the chromogenic solution contains TMB and H₂O₂. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18;980-985.
2. HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
3. ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
4. WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
5. WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
6. AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
7. BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
8. NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
9. PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
10. QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (µl)	SAMPLE(S) CONTROLS (µl)
Calibrators (0-5) Controls, Samples Working Anti-CT-HRP conjugate	100 - 50	- 100 50
Incubate for 18 ± 1 hours at 2 – 8°C. Aspirate the contents of each well. Wash 3 times with 400 µl of Wash Solution and aspirate.		
Chromogenic TMB Solution	100	100
Incubate for 30 min at room temperature		
Stop Solution	100	100
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (versus 630 or 650 nm)		

DIAsource Catalogue Nr : KAP0421	P.I. Number : 1700468/en	Revision nr : 110221/1
-------------------------------------	-----------------------------	---------------------------

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XVII. BIBLIOGRAPHIE

1. GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVISKINER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18;980-985.
2. HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
3. ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
4. WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
5. WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
6. AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
7. BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
8. NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
9. PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
10. QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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2. HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.

3. ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.

4. WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.

5. WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.

6. AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.

7. BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.

8. NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.

9. PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.

10. QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57

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XVII. LITERATUR

1. GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVISKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18;980-985.
2. HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974) **A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.**
J. Clin. Endocrinol. Metab., 39:487-495.
3. ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
4. WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
5. WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
6. AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
7. BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
8. NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
9. PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
10. QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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1. GRAZE K., SPILER I.J., TASHJIAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18;980-985.
2. HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
3. ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
4. WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
5. WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
6. AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and Foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
7. BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
8. NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
9. PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
10. QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.

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- GRAZE K., SPILER I.J., TASHUAN A.H., MELVIN K.E.W., CERVISKINER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978) **Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.** Engl. J. Med., 299,18;980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974) **A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.** J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983) **the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.** Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981) **Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.** J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)

Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.

Ann. Surg., 188,2:139-141.

- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985) **Parathyroid hormone, calcitonin, and the calciferols.** In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987) **SCC antigen and other tumor markers in lung cancer: preliminary results.** Excerpta Medica, 162-170.
- NICOLI P. et al. (1995) **Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.** Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994) **Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.** J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994) **Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.** J. Bone Miner Res. 9, 1, 53-57.

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
PREC AGENT	Precipitating Agent
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
EXTR SOLN CONC	Extraction Buffer Concentrate
CART	Cartridge
SAV HRP	Streptavidin HRP
PIPETTE	Pipette
WASH SOLN	Wash buffer

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