# **Androstenedione ELISA**

***KAPD3265***

---

**LOT** : 140710/1

# ANDROSTENEDIONE ELISA

en

KAPD3265  
IN VITRO DIAGNOSTIC USE

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

## 1 INTRODUCTION

The **DIAsource Androstenedione Enzyme Immunoassay Kit** provides materials for the quantitative determination of Androstenedione in serum and EDTA plasma.

**This assay is intended for in vitro diagnostic use only.**

The steroid hormone Androstenedione is one of the main androgens, besides Testosterone and Dehydroepiandrosterone. Testosterone, the most important biological active androgen, is derived from peripheral enzymatic conversion of Androstenedione.

In males, androgens are secreted primarily by the Leydig cells of the testes, to some degree also in the adrenal cortex. In females, the androgens are secreted mainly in the adrenal glands and in the ovary. Around 10% of the androgens are derived from peripheral conversion, mainly of DHEA. Androstenedione and Testosterone show high diurnal variability. The highest levels are measured in the morning. At the age of puberty serum androstenedione levels rise, after menopause they decline again. High androstenedione levels are measured during pregnancy. In women, high levels of androstenedione (47-100% above normal) are generally found in hirsutism, mostly in combination with other androgens as testosterone and DHEA-S. Androstenedione overproduction is due to ovarian dysfunction or maybe of adrenal origin. High circulating androstenedione levels are found in women with polycystic ovaries and 21-hydroxylase effect. Significant lower androstenedione levels are found in postmenopausal osteoporosis.

## 2 PRINCIPLE OF THE TEST

The DIAsource Androstenedione ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding.

The microtiter wells are coated with an antibody directed towards an antigenic site on the Androstenedione molecule. Endogenous Androstenedione of a patient sample competes with an Androstenedione horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase conjugate is reverse proportional to the concentration of Androstenedione in the sample. After addition of the substrate solution, the intensity of colour developed is reverse proportional to the concentration of Androstenedione in the patient sample.

## 3 PRECAUTIONS

- This kit is for in vitro diagnostic use only. For professional use only.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
- The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
- Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
- Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Allow the reagents to reach room temperature (21-26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Avoid contact with *Stop Solution* containing 0.5 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
- Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
- TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
- For information on hazardous substances included in the kit please refer to Material Safety Data Sheets. Material Safety Data Sheets for this product are available upon request directly from DIAsource.

## 4 KIT COMPONENTS

### 4.1 Contents of the Kit

- |    |
|----|
| LL |
|----|

**Microtiterwells**, 12x8 (break apart) strips, 96 wells  
Wells coated with a polyclonal anti-Androstenedione antibody
- |     |   |
|-----|---|
| CAL | N |
|-----|---|

**Calibrators N= 0 to 5**, 6 vials, 1 mL, ready to use  
See exact values on vial labels  
Conversion: ng/mL x 3.492 = nmol/l,  
Standards are calibrated against Reference Material of the National Measurement Institute Australia (NMIA-M955)  
contain non-mercury preservative.
- |         |   |
|---------|---|
| CONTROL | N |
|---------|---|

**Control Low & High**, 2 vials, 1.0 mL each, ready to use;  
For control values and ranges please refer to vial label or QC-Datasheet.  
Contain non-mercury preservative.
- |    |     |
|----|-----|
| Ag | HRP |
|----|-----|

**Enzyme Conjugate**, 1 vial, 25 mL, ready to use  
Androstenedione conjugated to horseradish Peroxidase,  
\* contains non-mercury preservative.
- |       |     |
|-------|-----|
| CHROM | TMB |
|-------|-----|

**Substrate Solution**, 1 vial, 25 mL, ready to use  
TMB
- |      |      |
|------|------|
| STOP | SOLN |
|------|------|

**Stop Solution**, 1 vial, 14 mL, ready to use  
contains 0.5M H<sub>2</sub>SO<sub>4</sub>  
Avoid contact with the stop solution. It may cause skin irritations and burns.
- |      |      |      |
|------|------|------|
| WASH | SOLN | CONC |
|------|------|------|

**Wash Solution**, 1 vial, 30 mL (40X concentrated)  
see „Preparation of Reagents“

### 4.2 Equipment and material required but not provided

- A microtiterplate calibrated reader (450±10 nm)
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Distilled or deionized water.
- Timer
- Semi-logarithmic graph paper or software for data reduction

### 4.3 Storage and stability of the Kit

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2-8°C. Microtiter wells must be stored at 2-8°C. Once the foil bag has been opened, care should be taken to close it tightly again. Opened kits retain activity for three months if stored as described above.

### 4.4 Preparation of Reagents

Allow all reagents and required number of strips to reach room temperature prior to use.

#### Wash Solution

Add deionized water to the 40X concentrated Wash Solution

Dilute 30 mL of concentrated Wash Solution with 1170 mL deionized water to a final volume of 1200 mL.

*The diluted Wash Solution is stable for 2 weeks at room temperature.*

### 4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheets.

### 4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DIAsource ImmunoAssays S.A. has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

## 5 SPECIMEN

Serum or EDTA plasma can be used in this assay.

Do not use Heparin or Citrate plasma. Heparin plasma leads to slightly reduced values. For citrate plasma the results are significantly increased. Do not use haemolytic, icteric or lipaemic specimens.

Please note : Samples containing sodium azide should not be used in the assay.

### 5.1 Specimen Collection

#### Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

#### Plasma:

Whole blood should be collected into centrifuge tubes containing anti coagulant and centrifuged immediately after collection. (E.g. Sarstedt Monovette for EDTA plasma)

### 5.2 Specimen Storage

Specimens should be capped and may be stored for up to 5 days at 2-8°C prior to assaying.

Specimens held for a longer time should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.

### 5.3 Specimen Dilution

If in an initial assay, a serum specimen is found to contain more than the highest calibrator, the specimens can be diluted 10-fold or 100 fold with *Calibrator 0* and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

#### Example:

- a) dilution 1:10: 10 µL Sample + 90 µL Calibrator 0 (mix thoroughly)
- b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL Calibrator 0 (mix thoroughly).

## 6 TEST PROCEDURE

### 6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each calibrator, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

### 6.2 Assay Procedure

Each run must include a calibration curve.

1. Secure the desired number of Microtiterwells in the holder.
2. Dispense **20 µL** of each Calibrators, controls and samples with new disposable tips into appropriate wells.
3. Dispense **200 µL** Enzyme Conjugate into each well.
4. Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
5. Incubate for **60 minutes** at room temperature.
6. Briskly shake out the contents of the wells.  
Rinse the wells 4 times with diluted Wash Solution (400 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.  
**Important note:**  
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
7. Add **200 µL** of Substrate Solution to each well.
8. Incubate for **30 minutes** at room temperature.
9. Stop the enzymatic reaction by adding **100 µL** of Stop Solution to each well.
10. Read the OD at **450±10 nm** with a microtiter plate reader **within 10 minutes** after adding the Stop Solution.

### 6.3 Calculation of Results

1. Calculate the average absorbance values for each set of calibrators, controls and patient samples.
2. Using semi-logarithmic graph paper, construct a calibration curve by plotting the mean absorbance obtained from each calibrator against its concentration with absorbance value on the vertical(Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the calibration curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this calibration curve. Samples with concentrations higher than that of the highest calibrator have to be further diluted. For the calculation of the concentrations this dilution factor has to be taken into account.

### 6.3.1 Example of Typical Calibration curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Calibrator	Optical Units (450 nm)
Calibrator 0 (0 ng/mL)	1.93
Calibrator 1 (0.1 ng/mL)	1.76
Calibrator 2 (0.3 ng/mL)	1.39
Calibrator 3 (1.0 ng/mL)	0.87
Calibrator 4 (3.0 ng/mL)	0.47
Calibrator 5 (10.0 ng/mL)	0.23

## 7 EXPECTED VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy adults, using the DIAsource Androstenedione ELISA the following values are observed:

### Females

Age (years)	n	Mean (ng/mL)	2.5 <sup>th</sup> Percentile (ng/mL)	97.5 <sup>th</sup> Percentile (ng/mL)
0-10	29	0.39	0.02	0.86
11-17	17	1.36	0.25	2.78
18-53	66	2.22	0.75	3.89
54-82	26	1.32	0.35	2.49

### Males

Age (years)	n	Mean (ng/mL)	2.5 <sup>th</sup> Percentile (ng/mL)	97.5 <sup>th</sup> Percentile (ng/mL)
0-10	34	0.40	0.01	1.31
11-17	16	1.75	0.33	3.30
18-53	36	2.15	0.45	4.20
54-82	44	1.95	0.30	3.93

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

## 8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DIAsource directly.

## 9 ASSAY CHARACTERISTICS

### 9.1 Assay Dynamic Range

The range of the assay is between 0.021 – 10 ng/mL.

### 9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross-reactivity of the assay:

Compound	Crossreactivity %
Androstenedione	100.0
Androsterone	< 0.01
Aldosterone	0.0
Cortisol	< 0.2
Dihydrotestosterone	< 0.01
Dihydroepiandrosterone	0.01
Estriol	1.8
16-Epiestriol	< 0.01
Estradiol	< 0.01
Estriol-3-glucuronide	< 0.01
Estriol-16-glucuronide	< 0.01
Estriol-16-sulfate	< 0.01
Estrone	< 0.01
17a-Pregnenolone	< 0.01
17OH-Progesterone	< 0.3
Progesterone	< 0.01
Testosterone	0.6

### 9.3 Analytical Sensitivity

The analytical sensitivity was calculated from the mean minus two standard deviations of twenty (20) replicate analyses of *Calibrator0* and was found to be 0.021 ng/mL.

### 9.4 Precision

The within assay variability (Intra Assay) and between assay variability (Inter Assay) are shown below:

#### 9.4.1 Intra Assay Variation

Sample	Intra Assay Variation		
	n	Mean (ng/mL)	CV (%)
1	20	1.2	5.3
2	20	0.7	9.2
3	20	8	6.9

#### 9.4.2 Inter Assay Variation

Sample	Inter Assay Variation		
	n	Mean (ng/mL)	CV (%)
1	40	4.3	8.3
2	40	3.0	8.1
3	40	1.0	8.7

#### 9.4.3 Inter-Lot

The inter-lot variation was determined by repeated measurements of 3 samples in 6 replicates per sample with 3 different kit lots:

Sample	n	Mean Lot 1 (ng/mL)	Mean Lot 2 (ng/mL)	Mean Lot 3 (ng/mL)	Inter Lot CV (%)
1	18	5.77	5.25	4.85	8.7
2	18	3.43	3.17	3.50	5.2
3	18	7.12	7.42	6.63	5.6

### 9.5 Recovery

Samples have been spiked by adding Androstenedione solutions with known concentrations in a 1:1 ratio. The expected values were calculated by addition of half of the values determined for the undiluted samples and half of the values of the known solutions. The % Recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100.

		Sample 1	Sample 2	Sample 3
Concentration [ng/mL]		1.47	1.43	1.95
Average Recovery		106.6	102.6	104.0
Range of Recovery [%]	from	96.4	99.1	93.3
	to	112.9	108.6	114.0

### 9.6 Linearity

		Sample 1	Sample 2	Sample 3
Concentration [ng/mL]		6.22	4.70	7.20
Average Recovery		97.0	96.3	93.4
Range of Recovery [%]	from	92.3	91.6	86.7
	to	102.9	101.3	105.6

## 10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

### 10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

### 10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of Androstenedione in a sample.

### 10.3 High-Dose-Hook Effect

No hook effect was observed in this test.

## 11 LEGAL ASPECTS

### 11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test. The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DIAsource.

### 11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutical consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutical consequences.

### 11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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Revision date : 2014-07-10

	<b>Used symbols</b>
	Consult instructions for use
	Storage temperature
	Use by
<b>LOT</b>	Batch code
<b>REF</b>	Catalogue number
<b>CONTROL</b>	Control
<b>I V D</b>	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
<b>U U</b>	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
<b>STRIP</b>	Dot Strip
SUB	Substrate
EXTR SOLN CONC	Extraction Buffer Concentrate
<b>CART</b>	Cartridge
SAV HRP	Streptavidin HRP
PIPETTE	Pipette
WASH SOLN	Wash buffer

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