



INTENDED USE

Anti-Phospholipid 10 Dot is used for the qualitative detection of IgG or IgM antibodies to phospholipids and serum proteins in human serum, for the diagnosis of anti-phospholipid antibody syndrome (APAS).

APAS is an autoimmune disorder comprising such clinical symptoms as arterial or venous thrombosis, thrombocytopenia and recurrent fetal loss. Primary APAS, as well as systemic lupus erythematosus (SLE), is characterized by the appearance of autoantibodies to negatively charged phospholipids, including cardiolipin antibodies (1). Although the significance and pathological relevance of phospholipid antibodies are not yet completely elucidated, the detection of such autoantibodies is widely established and plays an important role in the diagnosis of systemic autoimmune diseases.

Anti-cardiolipin antibodies are most frequently used for the diagnosis of APAS. Antibodies to the other phospholipids are especially of importance for differential diagnosis, if anti-cardiolipin antibodies are negative. Anti-phospholipid antibodies of APAS patients bind preferentially to a complex of phospholipids and the co-factor beta-2-glycoprotein I (beta-2-GP I, Apolipoprotein H) (2, 3), whereas antibodies directed against phospholipids only can also be detected in patients suffering from infectious diseases. beta-2-GP-I, a serum protein with a molecular weight of 50 kDa, affects platelet aggregation and coagulation. The positively charged fifth domain of beta-2 GP-I interacts with negatively charged phospholipids. This interaction results in conformational changes of the protein and the creation of new epitopes, apparently recognized by autoimmune phospholipid autoantibodies. Beside this, antibodies reacting directly with serum proteins, e.g. annexin V, prothrombin, can be detected in APAS patients.

- (1) Harris EN, Gharavi AE, Boey ML, Patel BM, Mackworth-Young GG, Loizou S and Hughes GRV: Anticardiolipin antibodies: detection by radioimmunoassay and association with thrombosis in systemic lupus erythematosus. Lancet 1983 11:1211
(2) Galli M, Comfurius P, Maassen C, Hemker HC, DeBaets MHVan Breda-Vriesman PJC, Barbui T, Zwaal RFA, Bevers EM: Anticardiolipin antibodies (ACA) directed not to cardiolipin but to a plasma protein factor. Lancet 1990 335:1544-1547
(3) McNeil HP, Simpson RJ, Chesterman CN, Krilis SA: Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding factor of coagulation: beta 2-glycoprotein I (apolipoprotein H). Proc Natl Acad Sci USA 1990 87:4120-4124

Anti-Phospholipid 10 Dot

- 20 x 10 determinations -

IVD In-vitro diagnostic device



Enzyme immunodot for the determination of IgG- or IgM-antibodies to phospholipids and serum proteins in human serum

Table with 2 columns: REF (Catalogue number) and LOT (Batch code). Rows include: Consult accompanying documents, Temperature limitations, Consult operating instruction, Manufactured by, Use by, Biological risk.

PRINCIPLE of the TEST

Anti-Phospholipid 10 Dot is an immunodot assay used for the qualitative detection of IgG or IgM antibodies to phospholipids and serum proteins in human serum.

Each kit includes 20 numbered strips. The strips consist of a membrane where 10 different autoantigen lines and a positive control line are sprayed on. Each line contains highly purified antigen preparation: (cardiolipin, phosphatidic acid, phosphatidyl-choline, -ethanolamine, -glycerole, -inositole, -serine, annexin V, beta2GPI, prothrombin)

During the first incubation autoantibodies of the patient sample bind to the target antigen immobilized on the solid phase (membrane). Unbound sample components are removed by a washing step after an incubation period of 30 minutes at room temperature (RT) while shaking.

During a second step bound antibodies react specifically with anti-human-IgG or anti-human-IgM antibodies which are conjugated to horseradish peroxidase (POD). Excessive conjugate is removed from the solid-phase immune complex by an additional washing step after an incubation period of 15 minutes at RT while shaking.

The horseradish peroxidase converts the colourless substrate solution into a dark purple precipitating line on the membrane. After 10 minutes at RT while shaking the reaction is stopped by two washing steps.

After the strips have been dried the results can be read by eye.



GA GENERIC ASSAYS GmbH

Ludwig-Erhard-Ring 3

15827 Dahlewitz, Germany

Telefon: +49 (0) 33708-9286-0
Fax: +49 (0) 33708-9286-50

www.genericassays.com

## PATIENT SAMPLES

### Specimen collection and storage

Blood is taken by venipuncture. Serum is separated after clotting by centrifugation.

The samples may be kept at 2...8 °C for up to three days. Long-term storage requires - 20 °C.

Repeated freezing and thawing should be avoided. If samples are to be used for several assays, initially aliquot samples and keep them at - 20 °C.

### Preparation before use

Samples to be measured with the Anti-Phospholipid 10 Dot assay are used undiluted. Take care to agitate samples gently before pipetting to ensure homogeneity.

## TEST COMPONENTS for 20 x 10 determinations

<b>A</b>	<b>Dot strips</b>	20 dot strips vacuum sealed
<b>Ag</b>	20 strips with 11 test dot lines 10 test lines coated with highly purified antigens: - cardiolipin - phosphatidic acid - phosphatidylcholine - phosphatidylethanolamine - phosphatidylglycerol - phosphatidylinositol - phosphatidylserine - annexin V - β2GPI - prothrombin positive control	
<b>B</b>	<b>Buffer, 10-fold</b>	<b>15 ml</b>
<b>BUF</b>	sufficient for 150 ml	concentrate capped white
	<b>10x</b>	
<b>C</b>	<b>IgG conjugate, 21-fold</b>	<b>1.2 ml</b>
<b>CONJ</b> <b>G</b>	anti-human-IgG (rabbit) coupled with horseradish peroxidase (POD)	ready for use capped red
<b>D</b>	<b>IgM conjugate, 21-fold</b>	<b>1.2 ml</b>
<b>CONJ</b> <b>M</b>	anti-human-IgM (rabbit) coupled with horseradish peroxidase (POD)	ready for use capped green
<b>E</b>	<b>Substrate</b>	<b>11 ml</b>
<b>SOLN</b> <b>TMB</b>	3,3',5,5'-Tetramethylbenzidine	ready for use capped blue
<b>F</b>	<b>Incubation tray for 12 dot strips</b>	2 x
<b>G</b>	<b>Lot specific interpretation template</b>	1 x

### Materials required in addition

- micropipette 100 - 1000 µl
- micropipette 10 - 100 µl
- pipette tips
- shaker (rocking shaker recommended)
- graduated cylinders
- distilled or de-ionized water
- plastic tweezers
- paper towel

### Size and storage

The Anti-Phospholipid 10 Dot has been designed for 20 determinations of IgG or IgM antibodies to phospholipids.

The expiry date of each component is reported on its respective label, that of the complete kit on the box label.

Upon receipt, all components of the Anti-Phospholipid 10 Dot have to be kept at 2 - 8 °C, preferably in the original kit box.

After opening of the vacuum sealing, the test dot strips are stable for 4 weeks when stored in the plastic foil bag at 2...8°C.

After opening all other kit reagents are stable for at least 4 weeks if properly stored at 2...8°C.

### Preparation before use

All wash steps are performed at room temperature (RT). Therefore, the reagents (buffer solution, substrate) must have RT in time.

The dot strips are vacuum sealed. A sufficient number of strips has to be cut off with a scalpel or a cutter from the retaining membrane. Unused dot strips have to be kept dry and stored in the additionally supplied plastic foil bag.

Dilute the 10 fold concentrated buffer with de-ionized or distilled water (1+9).

For each test strip 5 ml of buffer solution is required

#### Example:

15 ml concentrated buffer + 135 ml distilled water.

The prepared solution (diluted from B) is stable at 2 - 8 °C up to 30 days.

All other components are ready for use.

Avoid exposure of the substrate to light.

### Cleaning procedure of the incubation tray

After application incubate the incubation tray for 30 min with a detergent and rinse with water subsequently.

In the following step fill with any type of alcohol (methanol, propanol or ethanol), incubate on the rocking shaker for 30 min and subsequently rinse with water.

Clean the incubation tray with a cotton bud, rinse with water, and let it dry.

## ASSAY PROCEDURE

- Follow the instructions strictly and avoid any time shift.
- The whole assay has to be performed on a shaker (rocking shaker at 40 – 50 per minute is recommended)
- The entire assay runs at RT (18...25 °C). Ensure that the required reagents (conjugate, buffer solution, substrate) are at RT.

1. Take the reagents and a sufficient number of dot strips out of the box, mix the reagents gently.
2. Dispense **1 ml** of buffer solution (made of B) in each well.
3. Place the strips with the reactive side down into the respective wells.
4. Add **30 µl serum** (final dilution 1 + 33)
5. Incubate for **30 minutes** at RT while shaking
6. Decant (**Caution:** Turn the incubation tray over carefully and gently decant the buffer solution, any remaining liquid has to be removed with an absorbent paper). Wash **5 min at RT** with **1 ml** buffer solution (made of B) while shaking.
7. Pipette **1 ml** buffer solution (made of B) and add into the respective wells.  
IgG determination: **50 µl** conjugate C  
IgM determination: **50 µl** conjugate D
8. Incubate for **15 minutes** at RT while shaking.
9. Decant (**Caution:** Turn the incubation tray over carefully and gently decant the buffer solution, any remaining liquid has to be removed with an absorbent paper). Wash **5 min at RT** with **1 ml** buffer solution (made of B) while shaking.
10. Pipette **0.5 ml** substrate (E) into the respective wells.
11. Incubate for **10 minutes** at RT while shaking.
12. Decant and wash for **2 min at RT** with **1 ml** buffer solution (made of B) while shaking.
13. Decant (**Caution:** Turn the incubation tray over carefully and gently decant the buffer solution, any remaining liquid has to be removed with an absorbent paper). Wash **2 min at RT** with **1 ml distilled water** while shaking to stop the reaction.
14. Collect the dot strips from the wells and dry the membranes by pressing the reactive side of the strip briefly onto absorbent paper. After approximately **30 min** the strips can be interpreted.

## DATA PROCESSING


The evaluation of the test results is performed by means of the provided lot specific evaluation template. To do this the strips must have been dried and glued onto the template.

The **positive control line** must be positive in all cases. The coloration of the line ensures that the test has been run correctly and the kit components are not degraded. If the positive control line shows no coloration the results **cannot** be interpreted.

The test lines are coated with highly purified human antigens and detect specific antibody binding of the sample in the test.

The intensity of the bands given on the evaluation template serves as a cut-off of each single band for the positive - negative decision.

## REFERENCE VALUES



### Positive result:

A sample is considered to be **positive** in respect to one of the antibodies if the coloration of the test line shows a more intense coloration than the band on the evaluation template.

### Negative result:

A sample is considered to be **negative** in respect to one of the antibodies if the coloration of the test line shows the same or less intense coloration than the band on the evaluation template.

### Validation:

In order to interpret the results the test line of the positive control has to show a clear coloration.

## Limitations of Method

Healthy individuals should be tested negative by the Anti-Phospholipid 10 Dot. However, asymptomatic individuals can show a positive antibody reaction.

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis.

## CHARACTERISTIC ASSAY DATA

### Specificity and Sensitivity

The diagnostic specificity of the Anti-Phospholipid 10 Dot was obtained through the measurement of 50 normal sera. IgG detection showed 3 samples had at least one antibody (specificity 94%). IgM detections showed that 9 sera in total tested positive for at least one antibody (specificity 82%).

The diagnostic sensitivity of the Anti-Phospholipid 10 Dot was obtained through the measurement of 64 sera from APAS patients. IgG detections showed 43 sera (67.2%) with at least one antibody (individual sensitivities: cardiolipin 62.5%, phosphatidic acid 42.2%, P-choline 0, P-ethanolamine 0, P-glycerol 32.8%, P-inositol 45.3%, P-serin 59.4%, annexin V 4.7%, β2GPI 56.3%, prothrombin 14.1%). In IgM, antibodies were detected in 58 sera (90.6%) (individual sensitivities: cardiolipin 71.9%, phosphatidic acid 60.9%, P-choline 0, P-ethanolamine 57.8%, P-glycerol 34.4%, P-inositol 42.2%, P-serin 62.5%, annexin V 46.9%, β2GPI 65.36%, prothrombin 65.5%). In 61 sera, at least one IgG or IgM antibody was found (95.3%).

## INCUBATION – SCHEME

# Anti-Phospholipid 10 Dot (5012)

All steps are performed at RT (18...25°C); the required reagents (dot strips, buffer solution, conjugate) need to have RT in time.

1.	Bring all needed reagents and number of strips to room temperature (18...25 °C) before use	
2.	Pipette <b>1 ml</b> buffer solution (diluted from B) for each strip	
3.	Put strips with reaction side to the bottom on the buffer solution	
4.	Add <b>30 µl</b> patient serum afterwards: (final dilution 1+33)	
5.	Incubate	While shaking, <b>30 minutes</b> , RT
6.	Decant, washing of the strips	While shaking, RT, 5 minutes with 1 ml buffer solution (diluted from B)
7.	Pipette <b>1 ml</b> buffer solution (diluted from B) in each well of the tray, add: for IgG determination: <b>50 µl</b> conjugate C for IgM determination: <b>50 µl</b> conjugate D	
8.	Incubate	While shaking, <b>15 minutes</b> , RT
9.	Decant, washing of the strips	While shaking, RT, 5 minutes with 1 ml buffer solution (diluted from B)
10.	Pipette <b>0.5 ml</b> substrate (E) into the respective well	
11.	Incubate	While shaking, <b>10 minutes</b> , RT (18...25°C)
12.	Decant, washing of the strips	While shaking, RT, 2 minutes with <b>1 ml</b> buffer solution (diluted from B)
13.	Decant, washing of the strips, stopping of the reaction	While shaking, RT, 2 minutes with <b>1 ml</b> distilled water
14.	Dry test strips on paper towel for 30 minutes, read results	

## SAFETY PRECAUTIONS

- **This kit is for in vitro use only.** Follow the working instructions carefully. GA GENERIC ASSAYS GmbH and its authorized distributors shall not be liable for damages indirectly or consequentially brought about by changing or modifying the procedure indicated.
- The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for re-constituted reagents.
- Do not use or mix reagents from different lots.
- Do not use reagents from other manufacturers.
- Avoid time shift during pipetting of reagents.
- All reagents should be kept at 2 - 8 °C prior use in the original shipping container.
- Some of the reagents contain small amounts of kathon (1% v/v) as a preservative. They must not be swallowed or allowed to come into contact with skin or mucosa.
- Since the kit contains potentially hazardous materials the following precautions should be observed:
  - Do not smoke, eat or drink while handling kit material,
  - Always use protective gloves,
  - Never pipette material by mouth,
  - Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant.
- In any case GLP should be applied with all general and individual regulations to the use of this kit.