DHEA-S-ELISA

KAPD1562

LOT : 110225/1

DHEA-S ELISA

en

KAPD1562

IN VITRO DIAGNOSTIC USE

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

1 INTRODUCTION

The **DIAsource DHEA-S Enzyme Immunoassay Kit** provides materials for the quantitative determination of DHEA-S in serum and plasma. **This assay is intended for in vitro diagnostic use only.**

Dehydroepiandrosterone (5-Androstene-3 β -0L-17-one, Androstenolone, Dehydroisoandrosterone, Transdehydroandrosterone, DHEA) is a steroid hormone present in blood mostly in its sulfate form (DHEA-S).

DHEA-S is a more specific product of the adrenals and measurements of this steroid are widely used in clinical practice. The clinical importance of plasma assays of DHEA-S is associated with the diagnosis of adrenal hyperplasia and differential diagnosis of hirsutism (1).

2 PRINCIPLE OF THE TEST

The DIAsource DHEA-S ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. The microtiter wells are coated with a polyclonal antibody directed towards a unique antigenic site of the DHEA-S molecule.

Endogenous DHEA-S of a patient sample competes with a DHEA-S horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off.




The amount of bound peroxidase conjugate is reverse proportional to the concentration of DHEA-S in the sample. After addition of the substrate solution, the intensity of colour developed is reverse proportional to the concentration of DHEA-S in the patient sample.

3 PRECAUTIONS

- This kit is for in vitro diagnostic use only.
- For information on hazardous substances included in the kit please refer to Material Safety Data Sheets.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according the national biohazard safety guideline or regulation.

4 KIT COMPONENTS

4.1 Contents of the Kit

1.  12x8 (break apart) strips, 96 wells
Wells coated with anti-DHEA-S antibody
2.  CAL N
N=0 to 6
7 vials, 1 ml, ready to use
See exact value on the vial label
Conversion : 1 μ g/ml = 2.6 μ mol/l
Contains 0.03% Proclin 300 + 0.005% gentamicin sulfate as preservative
3.  Ag HRP
Enzyme Conjugate, 1 vial, 25 ml, ready to use
DHEA-S conjugated to horseradish peroxidase
*Contains 0.03% Proclin 300, 0.015% BND and 0.010% MIT as preservative
4.  CHROM TMB
Substrate Solution, 1 vial, 14 ml, ready to use
TMB
5.  STOP SOLN
Stop Solution, 1 vial, 14 ml, ready to use
contains 0.5M H₂SO₄
Avoid contact with the stop solution. It may cause skin irritations and burns.
6.  WASH SOLN CONC
Wash Solution, 1 vial, 30 ml (40X concentrated)
see „Preparation of Reagents“

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Note: Additional Zero Calibrator for Sample dilution available on request.

4.1.1 Equipment and material required but not provided

1. A microtiterplate calibrated reader (450±10 nm)
2. Calibrated variable precision micropipettes.
3. Absorbent paper.
4. Distilled water.

4.2 Storage and stability of the Kit

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. All opened reagents must be stored at 2-8°C. Microtiter wells must be stored at 2-8°C. Once the foil bag has been opened, care should be taken to close it tightly again.

4.3 Preparation of Reagents

Allow all reagents and required number of strips to reach room temperature prior to use.

Wash Solution

Add deionized water to the 40X concentrated Wash Solution

Dilute 30 ml of concentrated Wash Solution with 1170 ml deionized water to a final volume of 1200 ml. The diluted Wash Solution is stable for 2 weeks at room temperature.

4.4 Disposal of the Kit

The disposal of the kit must be made according to the national official regulations. Special information for this product are given in the Material Safety Data Sheets.

5 SPECIMEN

Serum or plasma (EDTA-, Heparin- or citrat plasma) can be used in this assay. Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

5.1 Specimen collection

Serum :

Collect blood by venipuncture (e.g Sarstedt Monovette # 02.1388.001), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma :

Whole blood should be collected into centrifuge tubes containing anti coagulant and centrifuged immediately after collection.

(E.g for EDTA plasma Sarstedt Monovette – red cap - # 02.166.001; for Heparin plasma Sarstedt Monovette – orange cap - # 02.165.001; for Citrat plasma Sarstedt Monovette – green cap - # 02.167.001.)

5.2 Specimen storage

Specimens should be capped and may be stored for up to 5 days at 2-8°C prior to assaying.

Specimens held for a longer time should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen dilution

If in an initial assay, a serum specimen is found to contain more than the highest calibrator, the specimens can be diluted with calibrator 0 and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10 : 10 µl Serum + 90 µl Calibrator 0 (mix thoroughly)
- b) dilution 1:100 : 10 µl dilution a) 1:10 + 90 µl Calibrator 0 (mix thoroughly).

6 TEST PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipet tips for each calibrator, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents be ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Assay Procedure

Each run must include a calibration curve.

1. Secure the desired number of Microtiterwells in the holder.
2. Dispense **25 µl** of each Calibrator, controls and samples with new disposable tips into appropriate wells.
3. Dispense **200 µl** Enzyme Conjugate into each well.
4. Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
5. Incubate for **60 minutes** at room temperature without covering the plate.
6. Briskly shake out the contents of the wells.
Rinse the wells 3 times with diluted Wash Solution (400 µl per well). Strike the wells sharply on absorbent paper to remove residual droplets.
Important note:
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
7. Add **100 µl** of Substrate Solution to each well.
8. Incubate for **15 minutes** at room temperature.
9. Stop the enzymatic reaction by adding **50 µl** of Stop Solution to each well.
10. Read the OD at **450±10 nm** with a microtiterplate reader **within 10 minutes** after adding the Stop Solution.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of calibrators, controls and patient samples.
2. Construct a calibration curve by plotting the mean absorbance obtained from each calibrator against its concentration with absorbance value on the vertical(Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the calibration curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this calibration curve. Samples with concentrations higher than that of the highest calibrator have to be further diluted. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Calibration Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Calibrator	Optical Units (450 nm)
Calibrator 0 (0 µg/ml)	1.69
Calibrator 1 (0.1 µg/ml)	1.35
Calibrator 2 (0.5 µg/ml)	0.93
Calibrator 3 (1.0 µg/ml)	0.67
Calibrator 4 (2.5 µg/ml)	0.46
Calibrator 5 (5.0 µg/ml)	0.33
Calibrator 6 (10.0 µg/ml)	0.23

7 EXPECTED VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy adults, using the DHEA-S ELISA the following values are observed:

Men	< 50 years	0.59 – 2.96 µg/mL
	> 50 years	mean: 1.01 µg/mL
Women	< 50 years	0.40 – 2.17 µg/mL
	> 50 years	mean: 0.63 µg/mL

8. QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DIAsource directly.

9. ASSAY CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0 – 10 µg/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

Steroid	% Crossreactivity
DHEA-S	100.0
Androstenedione	20.9
Androsterone	8.5
Androsterone-Sulfate	< 0.1
Progesterone	4.7
Testosterone	0.3
Estrone	0.9
Estriol	< 0.1
17-β-Estradiol	< 0.1
Estradiol Sulfate	< 0.1
Cortisol	< 0.1

9.3 Analytical Sensitivity

The analytical sensitivity was calculated from the mean minus two standard deviations of twenty (20) replicate analyses of *Calibrator 0* and was found to be 0.044 µg/mL.

9.4 Precision

9.4.1 Intra Assay Variation

The within assay variability is shown below:

Sample	n	Mean (µg/mL)	CV (%)
1	20	0.2	4.6
2	20	1.3	5.1
3	20	4.2	4.7

9.4.2 Inter Assay Variation

The between assay variability is shown below:

Sample	Mean (µg/mL)	CV (%)
1	0.2	6.5
2	1.3	10.6
3	4.4	5.4

9.5 Recovery

Samples have been spiked by adding DHEA-S solutions with known concentrations in a 1:1 ratio.

The expected values were calculated by addition of half of the values determined for the undiluted samples and half of the values of the known solutions. The % Recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100.

Sample	Added Concentration 1:1 (v/v) (µg/mL)	Measured Conc. (µg/mL)	Expected Conc. (µg/mL)	Recovery (%)
1	--	0.2	0.20	100
	1.0	0.7	0.60	109
	2.5	1.4	1.35	106
	5.0	2.8	2.60	108
2	--	1.1	1.10	100
	1.0	1.1	1.05	101
	2.5	2.1	1.80	114
	5.0	3.4	3.05	110
3	--	3.8	3.80	100
	2.5	2.8	3.15	91
	5.0	4.3	4.40	97
	10.0	7.0	6.90	102

9.6 Linearity

Sample	Dilution	Mean Conc. (µg/ml)	Recovery (%)
1	None	4.825	-
	1:2	2.505	104
	1:4	1.191	99
	1:8	0.644	107
	1:16	0.320	106
2	None	1.147	-
	1:2	0.620	108
	1:4	0.310	108
	1:8	0.153	107
	1:16	0.078	109
3	None	4.200	-
	1:2	2.116	101
	1:4	1.024	98
	1:8	0.486	93
	1:16	0.264	101

10. LIMITATIONS OF USE

Any improper handling of samples or modification of this test might influence the results.

Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of DHEA-S in a sample.

High-Dose-Hook Effect

No hook effect was observed in this test.

11. LEGAL ASPECTS

Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DIAsource.

Therapeutical Consequences

Therapeutical consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutical consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutical consequences.

Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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1. REFERENCIAS

1. Tietz, N. W., Textbook of Clinical Chemistry, Saunders, 1968

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µP	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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