IL-10-EASIA

KAP1321

LOT : 110221/3

Read entire protocol before use.

IL-10-EASIA

I. INTENDED USE

Immunoenzymetric assay for the in vitro quantitative measurement of human interleukin-10 (IL-10) in serum.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource IL-10-EASIA Kit
- B. **Catalogue number :** KAP1321 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. *Biological activities*

Human interleukin-10 (IL-10) is a 19 kDA lymphokine produced by T helper lymphocytes, by monocytes, macrophages and B-lymphocytes. IL-10 was first characterized as a cytokine synthesis inhibitory factor (CSIF) able to inhibit cytokine synthesis by TH1 clones activated in the presence of antigen presenting cells. However, in the absence of monocytes, IL-10 directly inhibits the growth of T-cells triggered by immobilized anti-CD3 MoAb. This proliferation inhibition was found to be a result of specific inhibition of IL-2 production by the responding T-cells. In vitro, IL-10 is a very powerful inhibitor of monokines (including TNF- α , IL-1, IL-6 and IL-8) produced by LPS-activated monocytes and macrophages. The addition of IL-10 to B lymphocytes results in limited cell proliferation but most importantly in very high immunoglobulin production, a result of the transformation of B-cells into plasma cells. Finally, natural killer (NK) cells appear to be another target for the anti-inflammatory properties of IL-10. Indeed, recent data have shown that IL-10 can inhibit antigen induced IFN- γ production by NK-cells by inhibiting not only production but also the stimulatory effects of IL-12 and TNF on IFN- γ production.









B. *Clinical application*

So far, circulating levels of IL-10 have been found in serum of patients suffering of Non-Hodgkin's lymphoma, multiple myeloma, cerebral malaria or septic shock.

IV. PRINCIPLES OF THE METHOD

The DIIASource IL-10-EASIA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplate. The assay uses monoclonal antibodies (MAbs) directed against distinct epitopes of IL-10. Calibrators and samples react with the capture monoclonal antibody (MAB 1) coated on microtiter well and with a monoclonal antibody (MAB 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAB 1 – human IL-10 – MAB 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the IL-10 concentration. A calibration curve is plotted and IL-10 concentration in samples is determined by interpolation from the calibration curve. The use of the EASIA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Color Code	Reconstitution
 Microtiterplate with 96 anti IL-10 (monoclonal antibodies) coated wells	96 wells	blue	Ready for use
 Conjugate: HRP labelled anti-IL-10 (monoclonal antibodies) in TRIS-Maleate buffer with bovine serum albumin and thymol	1 vial 6 ml	red	Ready for use
 Calibrator N = 0 to 5 (see exact values on vial labels) in human plasma with benzamidine and thymol	6 vials lyophil.	yellow	Add 1 ml distilled water
 Specimen Diluent: human plasma with benzamidine and thymol	3 vials lyophil.	black	Add distilled water (see on the label for the exact volume)
 Incubation Buffer: Phosphate buffer with bovine serum albumin and thymol	1 vial 11 ml	black	Ready for use
 Wash Solution (Tris-HCl)	1 vial 10 ml	brown	Dilute 200 x with distilled water (use a magnetic stirrer).
 Controls - N = 1 or 2 in human plasma with thymol	2 vials lyophil.	silver	Add 1 ml distilled water
 Chromogen TMB (Tetramethylbenzidine) in Dimethylformamide	1 vial 1 ml	green	Dilute 0.2 ml into 1 vial of substrate buffer
 Substrate buffer: H ₂ O ₂ in acetate / citrate buffer	3 vials 21 ml	white	Ready for use
 Stop Solution: H ₂ SO ₄ 1.8N	1 vial 6 ml	black	Ready for use

Note: 1. Use Specimen Diluent for sample dilutions.
2. 1 pg of the calibrator preparation is equivalent to 5 mIU of the NIBSC 1st RR 93/722.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 50 µl, 100 µl, 200 µl, 1 ml and 10 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Horizontal microtiterplate shaker capable of 700 rpm ± 100 rpm
- Washer for Microtiterplates
- Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (bichromatic reading)
- Optional equipment: The ELISA-AID™ necessary to read the plate according to polychromatic reading (see paragraph XIA.) can be purchased from Robert Maciels Associates, Inc. Mass. 0.2174 USA.

VII. REAGENT PREPARATION

- Calibrators:** Reconstitute calibrators with 1 ml distilled water.
- Controls:** Reconstitute the controls with 1 ml distilled water.
- Specimen Diluent:** Reconstitute specimen diluent to the volume specified on the vial label with distilled water
- Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.
- Revelation Solution:** pipette 0.2 ml of the chromogen TMB into one of the vials of substrate buffer (H₂O₂ in acetate/citrate buffer). Extemporaneous preparation is recommended.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, calibrators, controls and Specimen Diluent are stable for 4 days at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 2 months. Avoid successive freeze thaw cycles.
- The concentrated Wash Solution is stable at room temperature until expiration date.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, the conjugate is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- The freshly prepared revelation solution is stable, before use, for maximum 15 minutes at room temperature and must be discarded afterwards.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum must be removed as soon as possible from the clot of red cells after clotting and centrifugation, and kept at 4°C. If the samples are not used immediately, they must be kept at -20°C for maximum 2 months, and at -70°C for longer storage (maximum one year).
- Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.
- Sampling conditions can affect values, therefore, strict precautions have to be taken during sampling to avoid impurities contained in sampling materials that would stimulate IL-10 production by blood cells and thus falsely increase serum IL-10 values.
- Collection tubes must be pyrogen-free.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Revelation Solution and the Stop Solution avoid pipettes with metal parts.

High precision pipettes or automated pipetting equipment will improve the precision.

Respect the incubation times.

To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).

Prepare a calibration curve for each run, do not use data from previous runs.

The Revelation Solution should be colourless. If a blue colour develops within a few minutes after preparation, this indicates that the reagent is unusable, and must be discarded.

Dispense the Revelation Solution within 15 minutes following the washing of the microtiterplate.

During incubation with Revelation Solution, avoid direct sunlight on the microtiterplate.

B. Procedure

- Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
- Secure the strips into the holding frame.
- Pipette 100 µl of Incubation Buffer into all the wells
- Pipette 100 µl of each Calibrator, Control and Sample into the appropriate wells.
- Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.
- Aspirate the liquid from each well.
- Wash the plate 3 times by:
 - Dispensing 0.4 ml of Wash Solution into each well
 - Aspirating the content of each well
- Pipette 100 µl specimen diluent and then 50 µl of anti-IL-10-HRP conjugate into all the wells.
- Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.
- Aspirate the liquid from each well.
- Wash the plate 3 times by:
 - Dispensing 0.4 ml of Wash Solution into each well
 - Aspirating the content of each well
- Pipette 200 µl of the freshly prepared revelation solution into each well within 15 minutes following the washing step.
- Incubate the microtiterplate for 30 minutes at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm, avoid direct sunlight.
- Pipette 50 µl of Stop solution into each well.
- Read the absorbencies at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 3 hours and calculate the results as described in section XI.

XI. CALCULATION OF RESULTS

A. Polychromatic Reading:

- In this case, the ELISA-AID™ software will do the data processing.
- The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).
- A second reading is performed at 490 nm against the same reference filter.
- The ELISA-AID™ Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.
- The principle of polychromatic data processing is as follows:
 - $X_i = OD$ at 450 nm
 - $Y_i = OD$ at 490 nm
 - Using a standard unweighted linear regression, the parameters A & B are calculated: $Y = A \cdot X + B$
 - If $X_i < 3$ OD units, then X calculated = X_i
 - If $X_i > 3$ OD units, then X calculated = $(Y_i - B)/A$
 - A 4-parameter logistic curve fitting is used to build up the calibration curve.
 - The IL-10 concentration in samples is determined by interpolation on the calibration curve.

B. Bichromatic Reading

- Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
- Calculate the mean of duplicate determinations.
- On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of IL-10 (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
- Read the concentration for each control and sample by interpolation on the calibration curve.
- Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

IL-10-EASIA		OD units Polychromatic model
Calibrator	0 pg/ml	0.047
	20.5 pg/ml	0.119
	60.0 pg/ml	0.265
	204.0 pg/ml	0.726
	691.0 pg/ml	2.077
	1976.0 pg/ml	4.138

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 1.6 pg/ml.

B. Specificity

No significant cross-reaction was observed in presence of 50 ng of IL-1α, IL-1β, IL-1ra, IL-2, IL-3, IL-4, IL-6, IL-7, IL-8, TNF-α, TNF-β, IFN-α, IFN-β, IFN-γ, TGF-β, GM-CSF, OSM, MIP-1α, MIP-1β, LIF, MCP-1, G-CSF, RANTES, PF-4, βTG, GRO, IP-10 and SCF. This IL-10 assay is specific for human natural and recombinant IL-10.

A very low level (<0.2%) of cross-reaction was observed with BCRF1 (viral IL-10) at a concentration of 70000pg/ml. BCRF1 gave a signal corresponding to 134 pg/ml of IL-10.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> ± SD (pg/ml)	CV (%)	Serum	N	<X> ± SD (pg/ml)	CV (%)
Serum 1	24	86 ± 2.4	2.8	Serum 1	12	90 ± 2.5	2.8
2	24	324 ± 12	3.7	2	14	335 ± 9	2.7

SD: Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added IL-10 (pg/ml)	Recovered IL-10 (pg/ml)	Recovery (%)
serum	0	0	-
	60	56	93
	215	215	100
	760	745	98

DILUTION TEST

Sample	Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)
serum	1/1	910	910
	1/2	455	390
	1/4	228	213
	1/8	114	107
	1/16	57	57

Samples were diluted with Specimen Diluent.

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrators have been added to the coated wells.

TIME DELAY

sample	0 min	15 min	30 min
1	91	90	88
2	341	330	341
3	202	194	208
4	47	50	51
5	1141	1196	1228
6	284	294	297
7	136	133	137
8	263	263	291

F. Hook effect

A sample with an IL-10 concentration of 870 ng/ml gives a higher OD than the last calibrator point.

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls that contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

For guidance, the results of 32 serum samples from apparently healthy persons with low CRP levels, ranged between 0 and 3.3 pg/ml with a lean value of 0.2 pg/ml.

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains H₂SO₄, the chromogen contains TMB in Dimethylformamide, Substrate buffer contains H₂O₂. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. MOORE K.W., O'GARRA A., DE WAAL-MALEFYT R., VIEIRA P., MOSMANN T.R. (1993).
Interleukin-10
Annu. Rev. Immunol., 11:165-190.
2. BLAY J.Y., BURDIN N., ROUSSET F., LENOIR G., BIRON P., PHILIP T., BANCHEREAU J., FAVROT M.C. (1993).
Serum Interleukin-10 in Non-Hodgkin's Lymphoma : a prognostic factor.
Blood, 82:2169-2174.
3. BOGDAN C., BODOVOTZ Y., NATHAN C., (1991).
Macrophage deactivation by Interleukin-10.
J. Exp. Med., 174:1549-1555.
4. MELVILLE P., ROUSSET F., BANCHEREAU J., KLEIN B., BATAILLE R., (1992).
Serum Interleukin-10 in early stage multiple myeloma.
Lancet, 340:1544-1545.
5. PEYRON F., BURDIN N., RINGWALD P., VUILLEZ J.P., ROUSSET F., BANCHEREAU J., (1994).
High levels of circulating IL-10 in human malaria.
Clin. Exp. Immunol., 95:300-303.
6. MARCHANT A., DEVIERE J., BYL B., DE GROOTE D., VINCENT J.L., GOLDMAN M., (1994).
Interleukin-10 production during septicaemia.
Lancet, 343:707-708.
7. DE GROOTE D., MARCHANT A., FAUCHET F., JADOUL M., DEHART I., GERARD C., GEVAERT Y. et al. (1994).
Characterisation of monoclonal antibodies against human interleukin-10 and their use in an ELISA for the measurement of this cytokine.
J. of Immunol. Methods, 177:225-234.

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (µl)	SAMPLE(S) CONTROLS (µl)
Incubation Buffer	100	100
Calibrators (0-5)	100	-
Samples, Controls	-	100
Incubate for 2 hours at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 µl of Wash Solution and aspirate.		
Specimen Diluent	100	100
Anti-IL-10 -HRP conjugate	50	50
Incubate for 2 hours at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 µl of Wash Solution and aspirate.		
Revelation Solution	200	200
Incubate for 30 min at room temperature with continuous shaking at 700 rpm.		
Stop Solution	50	50
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (and 490 nm) versus 630 (or 650 nm)		

DIAsource Catalogue Nr : KAP1321	P.I. Number : 1702102/en	Revision nr : 110221/1
-------------------------------------	-----------------------------	---------------------------

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
PREC AGENT	Precipitating Agent
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
EXTR SOLN CONC	Extraction Buffer Concentrate
CART	Cartridge
SAV HRP	Streptavidin HRP
PIPETTE	Pipette
WASH SOLN	Wash buffer

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