Leptin-EASIA

KAP2281

LOT : 140606/2

Read entire protocol before use.

Leptin-EASIA

I. INTENDED USE

Immunoenzymetric assay for the in vitro quantitative measurement of human Leptin in serum.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource Leptin-EASIA Kit
- B. **Catalogue number :** KAP2281 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. Biological activities

Leptin, the product of the *ob* gene, is a hormone secreted by adipocytes (1). Animals with mutations in the *ob* gene are obese, diabetic and have reduced activity (2). Administration of recombinant leptin to these animals decreases food intake and causes weight loss. In humans, this type of mutation has not been found.

Human leptin cDNA encodes a 167 amino acid non-glycosylated protein including a 21 AA signal peptide, which is cleaved to give mature human leptin. The human receptor for leptin (OB-R) has been identified (7) as a 1144 amino acid transmembrane glycoprotein. It is expressed in the choroid plexus and in the hypothalamus. Leptin is implicated in an increasing number of endocrine regulations including adiposity (3), satiety, energy homeostasis (4), puberty and fertility (5, 6).

B. Clinical application






Little is known about the physiologic actions of leptin in humans. The availability of a sensitive and specific assay to measure concentrations of leptin in serum will accelerate our understanding.

IV. PRINCIPLES OF THE METHOD

The DIASource Leptin-EASIA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplate. The assay uses monoclonal antibodies (MAbs) directed against distinct epitopes of human Leptin. Calibrators and samples react with the capture monoclonal antibody (MAB 1) coated on microtiter well and with a monoclonal antibody (MAB 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAB 1 – human Leptin – MAB 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB ready for use) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the Leptin concentration.

A calibration curve is plotted and Leptin concentration in samples is determined by interpolation from the calibration curve. The use of the EASIA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Color Code	Reconstitution
 Microtiterplate with 96 anti Leptin (monoclonal antibodies) coated wells	96 wells	blue	Ready for use
 Conjugate: HRP labelled anti-Leptin (monoclonal antibodies) in TRIS-Maleate buffer with bovine serum albumin, thymol and ProClin 300 (< 0.05%)	1 vial 11 ml	red	Ready for use
 Zero Calibrator in bovine serum with thymol	1 vial lyophil.	yellow	Add 1 ml distilled water
 Calibrator N = 1 to 5 (see exact values on vial labels) in bovine serum with thymol	5 vials lyophil.	yellow	Add 0.5 ml distilled water
 Incubation Buffer: TRIS-Maleate buffer with bovine serum albumin, thymol and ProClin 300 (< 0.05%)	1 vial 6 ml	black	Ready for use
 Wash Solution (Tris-HCl)	1 vial 10 ml	brown	Dilute 200 x with distilled water (use a magnetic stirrer).
 Controls - N = 1 or 2 in human plasma with benzamidin and thymol	2 vials lyophil.	silver	Add 0.5 ml distilled water
 Chromogenic TMB Solution	1 vial 25 ml	white	Ready for use
 Stop Solution: HCl 2N	1 vial 25 ml	white	Ready for use

Note: 1. Use Zero Calibrator for sample dilutions.
2. 1 ng of the calibrator preparation is equivalent to 1 ng of the NIBSC 1st IS 97/594.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 50 µl, 100 µl, 500 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Horizontal microtiterplate shaker capable of 700 rpm ± 100 rpm

- Washer for Microtiterplates
- Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (bichromatic reading)
- Optional equipment: The ELISA-AID™ necessary to read the plate according to polychromatic reading (see paragraph XI.A.) can be purchased from Robert Maciels Associates, Inc. Mass. 0.2174 USA.

VII. REAGENT PREPARATION

- Calibrators:** Reconstitute the zero calibrator with 1 ml distilled water and the other calibrators with 0.5 ml distilled water.
- Controls:** Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, calibrators and controls should be frozen immediately after use and kept at -20°C for 3 months. Only one freeze thawing cycle is allowed, discard the calibrators and controls after the second use.
- The concentrated Wash Solution is stable at room temperature until expiration date.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, the conjugate is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum must be kept at 2 - 8°C.
- If the test is not run within 24 hours, storage in aliquots at -20°C is recommended. Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.
- Do not use haemolysed samples.
- Do not use lipemic samples.

X. PROCEDURE

- Handling notes**
Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.
High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).
Prepare a calibration curve for each run, do not use data from previous runs.
Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.
During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

B. Procedure

- Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
- Secure the strips into the holding frame.
- Pipette 50 µl of each Calibrator, Control and Sample into the appropriate wells.
- Pipette 100 µl of anti-Leptin-HRP conjugate into all the wells.
- Pipette 50 µl of Incubation Buffer into all the wells

6. Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm \pm 100 rpm.
7. Aspirate the liquid from each well.
8. Wash the plate 4 times by:
 - Dispensing 0.4 ml of Wash Solution into each well
 - Aspirating the content of each well
9. Pipette 100 μ l of Chromogenic Solution into each well within 15 minutes following the washing step.
10. Incubate the microtiterplate for 30 minutes at room temperature on a horizontal shaker set at 700 rpm \pm 100 rpm, avoid direct sunlight.
11. Pipette 200 μ l of Stop Solution into each well.
12. Read the absorbencies at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 3 hours and calculate the results as described in section XI.

XI. CALCULATION OF RESULTS

A. Polychromatic Reading:

1. In this case, the ELISA-AID™ software will do the data processing.
2. The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).
3. A second reading is performed at 490 nm against the same reference filter.
4. The ELISA-AID™ Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.
5. The principle of polychromatic data processing is as follows:
 - $X_i = OD$ at 450 nm
 - $Y_i = OD$ at 490 nm
 - Using a standard unweighted linear regression, the parameters A & B are calculated : $Y = A \cdot X + B$
 - If $X_i < 3$ OD units, then X calculated = X_i
 - If $X_i > 3$ OD units, then X calculated = $(Y_i - B)/A$
 - A 4-parameter logistic curve fitting is used to build up the calibration curve.
 - The Leptin concentration in samples is determined by interpolation on the calibration curve.

B. Bichromatic Reading

1. Read the plate at 405 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of Leptin (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
4. Read the concentration for each control and sample by interpolation on the calibration curve.
5. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

Leptin-EASIA		OD units
		Polychromatic model
Calibrator	0 ng/ml	0.019
	0.5 ng/ml	0.058
	1 ng/ml	0.106
	4.8 ng/ml	0.577
	27 ng/ml	3.052
	60 ng/ml	5.368

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 0.04 ng/ml.

B. Specificity

No significant cross-reaction was observed in presence of 3 μ g/ml of IL-1 α , IL-1 β , IL-4, IL-6, IL-8, IL-10, IL-15, TNF- α , TNF- β , IFN- γ , IP-10, ScF, GRO α , OSM, LIF, MCP-1, RANTES, PDGF, G-CSF, EGF, IL-12, PF4, MIG, C-peptide, IGF-1, IGF-2, Insulin & Glucagon.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> \pm SD (ng/ml)	CV (%)	Serum	N	<X> \pm SD (ng/ml)	CV (%)
A	10	1.5 \pm 0.2	13.3	A	20	5.9 \pm 0.6	10.2
B	10	9.0 \pm 0.9	10.0	B	20	18.9 \pm 2.4	12.7
C	10	43.4 \pm 1.5	3.5				

SD: Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added Leptin (ng/ml)	Theoretical Concent. (ng/ml)	Recovered Leptin (ng/ml)	Recovery (%)
Serum 1	0	7.2	7.6	105.6
	0.5	7.7	7.9	102.6
	2.4	9.6	9.1	94.8
	13.5	20.7	16.9	81.6
	30	37.2	34.5	92.7
Serum 2	0	13.3	12.4	93.2
	0.5	13.8	12.9	93.5
	2.4	15.7	14.5	92.4
	13.5	26.8	24.2	90.3
	30	43.3	38.6	89.1

DILUTION TEST

Sample	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
Serum 1	1/1	-	31.8
	1/2	15.9	15.6
	1/4	8.0	7.9
	1/8	4.0	3.9
	1/16	2.0	1.9
	1/32	1.0	0.8
Serum 2	1/1	-	32.8
	1/2	16.4	17.8
	1/4	8.2	8.3
	1/8	4.1	3.7
	1/16	2.1	1.7
	1/32	1.0	0.7

Samples were diluted with zero calibrator

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 40 minutes after the calibrators have been added to the coated wells.

TIME DELAY

Sample ng/ml	0 min	10 min	20 min	30 min	40 min
1	1.0	1.1	1.0	1.0	1.0
2	4.8	5.1	4.8	5.1	5.0
3	19.3	19.7	19.7	19.7	19.4
4	40.9	39.9	38.0	40.8	39.4

F. Hook effect

A sample spiked with 350 ng/ml Leptin gives higher OD's than the last calibrator point.

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.

- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls that contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Women (ages 16 to 80)				
BMI	n	Mean	SD	Range (ng/ml)
14-18	4	0.6	0.1	0.5 – 0.7
18-24	28	3.4	2.3	0.5 – 7.9
25-29	20	8.8	3.3	4.1 – 14.5
30-56	36	23.0	10.0	5.5 – 40.4

Men (ages 15 to 78)				
BMI	n	Mean	SD	Range (ng/ml)
14-18	1	0.5	-	-
18-24	20	1.2	0.9	0.5 – 3.2
25-29	27	4.2	3.9	0.5 – 14.6
30-56	23	10	10.5	2.5 – 42.1

The range is expressed as 2.5% to 97.5% percentiles.

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. FRIEMAN J. and al. (1998)
Leptin and regulation of body weight in mammals.
Nature ,395:763-769.
2. SMITH F.J. and al. (1998)
Brain administration of OB protein inhibits neuropeptide-Y- induced feeding of ob/ob mice.
Regulatory Peptides,75-76:433-439.
3. GUILLAUME M. and al. (1996)
Obesity in children,environmental and genetic aspects.
Horm. Metab. Res., 28 : 573-581.
4. KEIM N.L.and al. (1998)
Relation between circulating leptin cocentrations and appetite during a prolonged, moderate energy deficit in women.
American journal of Clinical Nutrition,68:794-801.
5. LEPERCQ J.and al. (1998)
Overexpression of placental leptin in diabetic pregnancy.
Diabetes,47:847-850.

6. MIYAKAWA M. and al.(1998)
Effect of growth hormone on serum concentrations of leptin : study in patients with acromegaly and GH deficiency.
The Journal of Clinical Endocrinology & Metabolism ,83:3476-3479.
7. TARTAGLIA L.A .and al.(1995)
Identification and expression cloning of a leptin receptor OR.
Cell , 83 (7):1263-1271.
8. CONSIDINE R..V. and al.(1996)
Serum immunoreactive-leptin concentrations in normal-weight and obese humans.
N. Engl. J. Med., 334(5):292-295.

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (µl)	SAMPLE(S) CONTROLS (µl)
Calibrators (0-5) Samples, Controls	50 -	- 50
Anti-Leptin-HRP Conjugate	100	100
Incubation Buffer	50	50
Incubate for 2 hours at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 4 times with 400 µl of Wash Solution and aspirate.		
Chromogenic Solution (TMB)	100	100
Incubate for 30 min at room temperature with continuous shaking at 700 rpm.		
Stop Solution	200	200
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (and 490 nm) versus 630 (or 650 nm)		

DIAsource Catalogue Nr : KAP2281	P.I. Number : 1700656/en	Revision nr : 140606/1
-------------------------------------	-----------------------------	---------------------------

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- FRIEMAN J. and al. (1998)
Leptin and regulation of body weight in mammals.
Nature, 395:763-769.
- SMITH F.J. and al. (1998)
Brain administration of OB protein inhibits neuropeptide-Y-induced feeding of ob/ob mice.
Regulatory Peptides, 75-76:433-439.

- GUILLAUME M. and al. (1996)
Obesity in children, environmental and genetic aspects.
Horm. Metab. Res., 28: 573-581.
- KEIM N.L. and al. (1998)
Relation between circulating leptin concentrations and appetite during a prolonged, moderate energy deficit in women.
American journal of Clinical Nutrition, 68:794-801.
- LEPERCQ J. and al. (1998)
Overexpression of placental leptin in diabetic pregnancy.
Diabetes, 47:847-850.
- MIYAKAWA M. and al. (1998)
Effect of growth hormone on serum concentrations of leptin: study in patients with acromegaly and GH deficiency.
The Journal of Clinical Endocrinology & Metabolism, 83:3476-3479.
- TARTAGLIA L.A. and al. (1995)
Identification and expression cloning of a leptin receptor OR.
Cell, 83 (7):1263-1271.
- CONSIDINE R.V. and al. (1996)
Serum immunoreactive-leptin concentrations in normal-weight and obese humans.
N. Engl. J. Med., 334(5):292-295.

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