

## Autoantibodies to intrinsic factor

**Intrinsic factor (IF)**, a sialic acid containing 60 kD glycoprotein, plays a vital role in the transport and absorption of vitamin B<sub>12</sub> in the intestine. After secretion by parietal cells of the stomach mucosa IF binds to **vitamin B<sub>12</sub>** ingested with food. This vitamin B<sub>12</sub>-IF complex is carried to the intestine where it allows absorption of vitamin B<sub>12</sub> via binding to a **specific IF receptor**. Subsequently vitamin B<sub>12</sub> is released into the blood binding to another protein (transcobalmin).

Reduced production of the IF and / or impairment of its transport function in the digestive system bring about a deficiency in vitamin B<sub>12</sub> leading to the development of **pernicious anemia** (Biermer's anemia). Patients suffering from atrophic **chronic gastritis of type A** exhibit autoantibodies to both parietal cell H<sup>+</sup>/K<sup>+</sup>-ATPase and the intrinsic factor produced by parietal cells. According to their binding sites these autoantibodies are divided into two types.

**Type 1 antibody** (blocking antibody) prevents the attachment of vitamin B<sub>12</sub> to IF in the stomach. In contrary **type 2 antibody** (binding antibody) interacts with both IF and IF-vitamin B<sub>12</sub> complex and prevents their absorption in the intestine by reacting with the IF region recognized by the mucosa.

### Etiopathogenesis of pernicious anemia

Pernicious anemia is the result of chronic atrophic gastritis of type A. Contrary to gastritis type B this disease is an autoimmune process with a progressive destruction of the gastric mucosa. Whilst type A gastritis involves the fundus and corpus of the stomach, type B gastritis affects the antrum as well and is usually associated with *Helicobacter pylori* infections.

#### Neurological complications

- Peripheral neuropathy
- Demyelination, axonal degeneration and neuronal death
- Sensory ataxia
- Memory loss
- Psychosis

### Clinical pattern

#### Histopathology

- Loss of gastric mucosa folds
- Submucosal infiltration of mononuclear cells
- Loss of parietal cells and replacement by mucus-containing cells
- Intestinal metaplasia

#### Progression

Typically type A gastritis progresses to an atrophic gastritis combined with pernicious anemia within 20-30 years.

Clinical symptoms of pernicious anemia occur in most cases after the age of 40.

#### Intestinal complications

- Patients with pernicious anemia face a higher risk to develop an adenocarcinoma.

#### Type A gastritis

- autoimmune pathogenesis
- fundus, corpus
- pernicious anemia
- autoantibodies to parietal cells (H<sup>+</sup>/K<sup>+</sup>-ATPase)
- autoantibodies to intrinsic factor
- achlorhydria
- low serum pepsinogen I
- hypergastrinemia

#### Type B gastritis

- non-autoimmune pathogenesis
- fundus, antrum, corpus
- *Helicobacter pylori* infection
- hypogastrinemia

## Diagnostic methods

- *Gastric biopsy*
- *Hematology*
  - Megaloblastic anemia
  - Macrocytosis with polymorphonuclear leucocytosis
- *Serology*
  - Low vitamin B<sub>12</sub> level
  - Antibodies to parietal cells
  - Antibodies to H<sup>+</sup>/K<sup>+</sup>-ATPase
  - Antibodies to intrinsic factor **(highly specific)**
  - Hypergastrinemia
  - Low level of serum pepsinogen I
- *Schilling test*
  - Confirmation of vitamin B<sub>12</sub> deficiency due to intrinsic factor deficiency

## Epidemiology

The overall prevalence of pernicious anemia has been estimated at 0.1% of the population. However, the incidence of pernicious anemia is age related.


A recent population study revealed that 1.9 percent of persons older than 60 years suffer from undiagnosed pernicious anemia (Carmel et al., 1996).

## Therapy

- Monthly intramuscular injection of 100 µg vitamin B<sub>12</sub>
- Daily oral intake of vitamin B<sub>12</sub> for people older than 60 years with diagnosed gastric atrophy

The detection of megaloblastic anemia and vitamin B<sub>12</sub> deficiency beside antibodies to intrinsic factor meet the criteria for differential diagnosis of pernicious anemia. Additional parameters (e.g. Schilling test) increase costs and expenditure of human labor.

**Detection of antibodies to parietal cells is not specific for pernicious anemia exclusively, because they are also detectable in other autoimmune diseases. Antibodies to the intrinsic factor are highly specific for pernicious anemia.**





# PERNICIOUS ANEMIA

## Anti-intrinsic factor (Order Code: 3600)

### Technical data:

<b>Assay</b>	96 determinations
<b>Disorder</b>	Pernicious anemia
<b>Parameters</b>	IgG antibodies to intrinsic factor (type 1 and type 2)
<b>Antigen</b>	Intrinsic factor (human, recombinant)
<b>Principle of the test</b>	Enzyme-linked immunosorbent assay (ELISA), semi-quantitative
<b>Results</b>	Binding index (BI), positive: $B > 1$ ; negative $BI < 0,8$
<b>Incubation scheme</b>	30 min RT – 30 min RT – 15 min RT
<b>Substrate</b>	Tetramethylbenzidine (TMB), 450 nm
<b>Sample</b>	human serum, 1 + 100

### Sensitivity and specificity

A specificity of 99% and a sensitivity of 98% has been established for GA Generic Assays GmbH Anti-intrinsic factor by investigating serum samples from healthy persons (n=110) and patients suffering from pernicious anemia (n=49)

### Precision

The coefficient of variation in 6 serums (15 pilot batches) was:

Intra-assay: < 7%  
Inter-assay: < 10%

### Literature

Carmel R. Prevalence of undiagnosed anemia in the elderly. Arch Intern Med 1996;156:1097-1100

Carmel R. Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anemia: influence of patient age and race. J Clin Exp Immunol 1992;89:74-77

Humbel RL, Pierrard V. Detection of anti-intrinsic factor antibodies by the ELISA technique. Aerztl Lab 1990;36:55-58

Tok BH, van Driel IR, Gleeson PA. Pernicious Anemia. N Engl J Med 1996;20:1441-1448

Waters HM, Dawson DW, Howarth JE, Geary CG. High incidence of type II autoantibodies in pernicious anaemia. J Clin Pathol 1993;46:45-57



# PERNICIOUS ANEMIA

Additional products for the differential diagnosis  
of pernicious anemia

Order Code	Product	Results	Determinations
3600	Anti-intrinsic factor (IF)	semi-quantitative	96
4020	BiermAK Dot (PCA, IF)	qualitative	24 x 2

- GA Generic Assays -

Your reliable partner in  
laboratory diagnostics!

Please contact for further information  
your local dealer!

Searching for more interesting products?

Please visit our homepage at:  
[www.genericassays.com](http://www.genericassays.com)