



Total Solution for DENGUE Test


Dengue COMBO Test, Dengue NS1 Antigen Test and Dengue IgG/IgM Test

Dengue Combo Test

Simultaneous Dengue NS1 + IgG/IgM Test

Excellent diagnostic tool for Dengue infection during all stage

- Complementary Dengue NS1 Ag, IgG and IgM
- No limitation of single test
- Fast result and Room temperature storage
- Available whole blood, plasma and serum
- High accuracy: 98.6% (NS1 Ag) and 98.6% (IgG/IgM)



Dengue IgG/IgM Ab Test

Dengue NS1 Ag Test



Dengue COMBO Test



CAT2123_Rev_00


For more information on our products, visit our website, www.humasis.com

Fast/Easy to use/Customer-oriented/Competitively-priced

For professional use

How to test

NS1 Ag Test




IgG/IgM Ab Test




Combo Test

NS1 Ag



IgG/IgM



Specifications

Product	CAT No.	Packing	Specimen	Result within	Storage	Shelf Life
Dengue Combo Test	ADEC-5025	25 Tests	W/B, Serum, Plasma NS1 Ag: 100µL IgG/IgM: 10µL	15~20 minutes	1~30°C (34~86°F)	24 Months
Dengue IgG/IgM test	ADEN-7025	25 Tests	100µL of W/B plasma/serum	15~20 minutes	1~30°C (34~86°F)	24 Months
Dengue NS1 Antigen Test	AEDG-7025	25 Tests	100µL of W/B plasma/serum	15~20 minutes	1~30°C (34~86°F)	24 Months

Ordering Information

Humasis Co., Ltd.
Rm.504, SHINWON VISION TOWER,
88, Jeonpa-ro, Dongan-gu, Anyang-si,
Gyeonggi-do, KOREA
TEL: +82-31-478-8591
FAX: +82-31-478-8586
Email: question@humasis.com
www.humasis.com

EC REP MT Promedt Consulting GmbH
Altenhofstr. 80
D-66386 St. Ingbert / Germany
TEL: +49 6894 - 58 10 20
FAX: +49 6894 - 58 10 21
Email: info@mt-procons.com
www.mt-procons.com

For more information on our products, visit our website, www.humasis.com