

# Disposabe Virus Sampling Tube

(Sterile)



**DAKEWE 达科为<sup>®</sup>**

**DAKEWE<sup>®</sup> | BioSci<sup>®</sup>**  
Shenzhen Dakewe Bio-engineering Co., Ltd.  
Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen  
Telephone: +86-0755-86235300 Website: <http://www.dakewe.com/>

OPEN

NMPA Reg. No.: Y5XB20200293

**Disposable Collection Swab** **STERILE**

- 1
- 2
- 3
- 4



## Kit component:

### # 3980212 MB

- Tube containing modi. viral transport medium 3ml
- Swab (Nasopharyngeal)

### # 3980222 MY

- Tube containing modi. viral transport medium 3ml
- Swab (Oropharyngeal )

# 3980212	MB
# 3980222	MY

Oropharyngeal (MY)



Nasopharyngeal (MB)



Containing Guanidine salt

Virus **inviabile**

Nucleotide purification required

PCR test

Anti bacteria and fungi

Sample stable up to 72 h at 2-8 °C

Stable at -70 °C



