

# Mole RT-qPCR PRODUCT LIST

info@molechina.com

NO	CAT.NO	PRODUCT DESCRIPTION	Specimen	Kit Size	Storage Temperature	Shelf Life	English IFU	Certificate
<b>Hepatitis</b>								
1	H012T048B0C1	Hepatitis B virus nucleic acid test kit (Quantitative)	Serum/Plasma	48 tests/kit	– (20±5) °C	12 months	Available	NMPA/ RUO
2	H051T025B0C0	Hepatitis A Virus Nucleic Acid Test Kit	Feces/Serum	25tests/kit	– (20±5) °C	9 months	Available	RUO
<b>Enterovirus</b>								
3	E011T048B0C1	EV71, CA16 and other enteroviral nucleic acid test kit	Feces/Anal swab/Pharynx swab	48 tests/kit	– (20±5) °C	9 months	Available	NMPA/ CE IVD RUO
4	E014T024B0C1	Enterovirus 71 (EV71) nucleic acid test kit	Feces/Anal swab	24 tests/kit	– (20±5) °C	-	N/A	NMPA/ RUO
5	E015T024B0C1	Coxsackie virus A16 (CA16) nucleic acid test kit	feces and oropharyngeal swab	24 tests/kit	– (20±5) °C	6 months	Available	NMPA/ RUO
6	E016T024B0C1	Enterovirus nucleic acid test kit	Feces, cerebrospinal fluid and oropharyngeal swab	24 tests/kit	– (20±5) °C	6 months	Available	NMPA/ RUO
7	E017T024B0C1	EV71 and CA16 Dual-channel Enterovirus test kit	Feces/Anal swab	24 tests/kit	– (20±5) °C	-	N/A	NMPA/ RUO
8	E061T050B0C0	Norovirus GI GII nucleic acid test kit	feces or vomit	50 tests/kit	– (20±5) °C	12 months	Available	RUO
9	E111T050B0C0	Coxsackie virus A6 (CA6) nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	Available	RUO
10	E112T050B0C0	Coxsackie virus A10 (CA10) nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	Available	RUO
11	E113T050B0C0	Coxsackie virus A6 and CA10 nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	Available	RUO
12	E121T050B0C0	Poliovirus type I nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	Available	RUO
13	E123T050B0C0	Poliovirus type II nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	Available	RUO
14	E124T050B0C0	Poliovirus type III nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	Available	RUO
15	E131T050B0C0	Sapovirus nucleic acid test kit	Feces or vomit	50 tests/kit	– (20±5) °C	9 months	Available	RUO
16	E141T050B0C0	Parachovirus nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	N/A	RUO
17	E151T050B0C0	Sapovirus, Astrovirus and Adenovirus nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	-	N/A	RUO
18	E152T050B0C0	Astrovirus nucleic acid test kit	feces or vomit	50 tests/kit	– (20±5) °C	12 months	Available	RUO
<b>Respiratory</b>								
19	R051T016B0C1	Chlamydia trachomatis, Chlamydia pneumoniae and Mycoplasma pneumoniae nucleic acid test kit	Sputum	16tests/kit	– (20±5) °C	12 months	Available	NMPA/ RUO
20	R052T024B0C1	Chlamydia Trachomatis Nucleic Acid Detection Kit	Sputum/Oropharyngeal swab/Urine	24 tests/kit	– (20±5) °C	12 months	Available	RUO
21	R061T016B0C1	<b>Mycoplasma pneumoniae and Macrolides-resistant strain nucleic acid test kit (A2063G, A2064G, P1 gene)</b>	★ Sputum	16tests/kit	– (20±5) °C	12 months	Available	NMPA/ CE IVD/ UKKA
22	R062T016B0C1	<b>Mycoplasma pneumoniae and Macrolides-resistant strain nucleic acid test kit (A2067G, C2617G)</b>	Sputum	16tests/kit	– (20±5) °C	12 months	Available	RUO
23	R070T050B0C0	Rubella Virus Nucleic Acid Quantitative Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
24	R071T050B0C0	Rubella virus nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
25	R072T050B0C0	Measles Virus nucleic acid test kit	Nasopharyngeal swab/Oropharyngeal swab/Sputum	50 tests/kit	– (20±5) °C	12 months	Available	RUO
26	R07301T2048	Rubella And Measles Nucleic Acid Test Kit	Oropharyngeal swabs, sputum and bronchoalveolar lavage fluid	48 Tests/Box	2-8 °C	24 months	Available	CE IVD
27	R075T050B0C0	Vaccine and wild measles virus typing kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
28	R081T050B0C0	Influenza A Virus Nucleic Acid Detection Kit	Oropharyngeal swab	50 tests/kit	– (20±5) °C	-	N/A	RUO
29	R082T050B0C0	Influenza B Virus Nucleic Acid Detection Kit	Oropharyngeal swab	50 tests/kit	– (20±5) °C	-	N/A	RUO
30	R083T050B0C1	<b>Influenza A and B virus nucleic acid test kit</b>	Oropharyngeal swab	50 tests/kit	– (20±5) °C	12 months	Available	NMPA/ CE IVD
31	R085T050B0C0	Influenza C virus nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
32	R091T050B0C0	Influenza A virus (seasonal H1) nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
33	R092T050B0C0	Influenza A virus (seasonal H3) nucleic acid test kit	Nasopharyngeal swab/Oropharyngeal swab	50 tests/kit	– (20±5) °C	9 months	Available	RUO
34	R093T050B0C0	Influenza A virus (seasonal H1 and H3) nucleic acid test kit	Nasopharyngeal swab/Oropharyngeal swab/Nasal swab	50 tests/kit	– (20±5) °C	9 months	Available	RUO
35	R101T050B0C0	Influenza A H1 Virus Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
36	R103T050B0C0	Influenza A virus (H1N1) nucleic acid test kit	Nasopharyngeal swab/Oropharyngeal swab	50 tests/kit	– (20±5) °C	9 months	Available	RUO
37	R111T050B0C0	Influenza B Yamagata Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
38	R112T050B0C0	Influenza B Victoria Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
39	R113T050B0C0	Influenza B Yamagata and Victoria test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
40	R121T050B0C0	Avian Influenza Virus H5 Subtype test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
41	R122T050B0C0	Influenza A Virus (H5N1) nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
42	R123T050B0C0	Avian Influenza Virus H7 Subtype test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
43	R124T050B0C0	Influenza A Virus (H7N9) nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
44	R125T050B0C0	Avian Influenza Virus H9 Subtype test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
45	R126T050B0C0	Influenza A Virus H5, H7 and H9 typing kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
46	R127T050B0C0	Avian Influenza Virus H4 Subtype test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
47	R129T050B0C0	Influenza A/H1N1 and A/H3, Influenza B/Yamagata and B/Victoria typing kit	nasopharyngeal swabs/oropharyngeal swabs/sputum	50 tests/kit	– (20±5) °C	9 months	Available	RUO
48	R131T032B0C1	Rhinovirus Nucleic Acid Detection Kit		32tests/kit	– (20±5) °C	-	N/A	RUO
49	R141T050B0C1	<b>EB Virus Nucleic Acid Detection Kit</b>	Whole Blood	50 tests/kit	– (20±5) °C	18 months	Available	NMPA/ RUO
50	R142T050B0C1	Respiratory Adenovirus Nucleic Acid Detection Kit	Sputum/Oropharyngeal swab	32/50 tests/kit	– (20±5) °C	12 months	Available	RUO

NO	CAT.NO	PRODUCT DESCRIPTION	Specimen	Kit Size	Storage Temperature	Shelf Life	English IFU	Certificate
51	R143T050B0C0	Adenoviruses, EB Viruses and cytomegaloviruses Nucleic Acid Multiplex test Panel	Sputum/ Sliva/ Nasopharyngeal swab/ Oropharyngeal swab /Bronchoalveolar lavage fluid	50 tests/kit	– (20±5) °C	12 months	Available	RUO
52	R144T050B0C0	Adenovirus types 3 and 7 typing kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
53	R145T050B0C0	Respiratory Adenovirus nucleic acid test kit (Quantitative)		50 tests/kit	– (20±5) °C	-	N/A	RUO
54	R146T050B0C0	Cytomegalovirus Nucleic Acid Test Kit	Respiratory swab, bronchoalveolar lavage fluid, saliva, and sputum	50 tests/kit	– (20±5) °C	12 months	Available	RUO
55	R151T050B0C1	Human Parainfluenza virus Nucleic Acid Detection Kit	nasopharyngeal swabs/oropharyngeal swabs/sputum	50 tests/kit	– (20±5) °C	6 months	Available	RUO
56	R152T050B0C0	Respiratory Parainfluenza Virus Type IV Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
57	R171T050B0C1	Bocavirus Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
58	R191T050B0C0	Middle East Respiratory Syndrome Nucleic Acid Detection Kit	nasopharyngeal swabs/oropharyngeal swabs/sputum	50 tests/kit	– (20±5) °C	6 months	Available	RUO
59	R201T032B0C0	Coronavirus SARS Nucleic Acid Detection Kit		32tests/kit	– (20±5) °C	-	N/A	RUO
60	R203T050B0C0	Coronavirus HKU1 Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
61	R205T002B0C2	2019-Novel Coronavirus (ORF1ab, E & N genes) (Reaction time: 70 mins, EXP:9 months, LOD: 500 copies/mL)			– (20±5) °C	9 months	Available	NMPA/CE IVD
62	R205T003B0C2 (last version)	2019-Novel Coronavirus (ORF1ab, E & N genes) (Reaction time: 90 mins, EXP:9 months, LOD: 500 copies/mL)			– (20±5) °C	9 months	Available	NMPA/CE IVD
63	R205T011B0C2 (new version)	2019-Novel Coronavirus (ORF1ab, E & N genes) (Reaction time: 90 mins, EXP: 12 months, LOD: 500 copies/mL)	Sputum/ Sliva/ Nasopharyngeal swab/ Oropharyngeal swab /Bronchoalveolar lavage fluid/ Feces	50 tests/kit 100 tests/kit 200 tests/kit	– (20±5) °C	12 months	Available	NMPA/CE IVD
64	R205T010B0C2	2019-Novel Coronavirus (ORF1ab, E & N genes) (Reaction time: 45 mins, only compatible with CFX-96, EXP:9 months, LOD: 100 copies/mL)			– (20±5) °C	9 months	Available	NMPA/CE IVD
65	R205T008B0C2	2019-Novel Coronavirus (ORF1ab, E & N genes) (Reaction time: 90 mins, Lyophilization, EXP: 24 months, LOD: 500 copies/mL)			<b>2-8°C (Lyophilized)</b>	<b>24 months</b>	Available	NMPA/CE IVD
66	R205T004B0C2	SARS-COV-2 Realtime PCR Test Kit (ORF1ab/N dualplex) (Reaction time: 70 mins, EXP:9 months, LOD: 500 copies/mL)	Sputum/ Sliva/ Nasopharyngeal swab/ Oropharyngeal swab /Bronchoalveolar lavage fluid/ Feces	50 tests/kit	– (20±5) °C	9 months	Available	NMPA/CE IVD
67	R208T050B0C2	2019-Novel Coronavirus, Influenza A and B virus nucleic acid test kit	Sputum/ Sliva/ Nasopharyngeal swab/ Oropharyngeal swab /Bronchoalveolar lavage fluid	50 tests/kit	– (20±5) °C	12 months	Available	CE IVD
68	R20802T2048	2019-Novel Coronavirus, Influenza A and B virus nucleic acid test kit	Oropharyngeal swab	48 tests/kit	<b>2-8°C (Lyophilized)</b>	<b>24 months</b>	Available	CE IVD
69	R213T050B0C2	2019-Novel Coronavirus, Influenza and RSV test kit	Sputum/Oropharyngeal swab/Bronchoalveolar lavage fluid	50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
70	R23102T2024	17 Multi Respiratory Pathogens test kit SARS-CoV-2 ORF1ab, Influenza A virus (IFA), Influenza B virus (IFB), Enterovirus (EV), human Metapneumovirus (hMPV A/B), Parainfluenza virus (PIV 1-4), Coronavirus (229E/OC43/NL63/HKU1), Respiratory Syncytial Virus (RSV A/B), Rhinovirus (RHV A/B/C), Adenovirus (ADV), Adenovirus 3, Adenovirus 7, EB virus (EB), Bocavirus (BoV), Bordetella pertussis, Mycoplasma Pneumoniae (MP), Legionella Pneumophila (LP)	Oropharyngeal swab	24 tests/kit	2-8°C (Lyophilized)	24 months	Available	CE IVD/ UKCA
71	R231T032B0C0	Human Respiratory Virus Detection Kit (18 types respiratory pathogens: Rhinovirus (RHV), respiratory syncytial virus (RSV), human metapneumovirus (hMPV), influenza A virus (IFA), influenza B virus (IFB), enterovirus (EV), parainfluenza virus 1-3 (PIV 1-3), bocavirus (BoV), cytomegalovirus (CMV), EB virus (EB), measles virus (MV), adenovirus (ADV) and coronavirus NL63/229E/OC43 and SARS	Oropharyngeal swab/Sputum	32 tests/kit	– (20±5) °C	9 months	Available	CE IVD
72	R23202T2012	Respiratory Microbes Nucleic Acid Test Kit (40 pathogens + 4 antibiotic resistant mutants)	swabs, sputum and alveolar lavage fluid	12 Tests/Box	2-8°C (Lyophilized)	12 months	Available	RUO
73	R241T050B0C0	Varicella-Zoster Virus Nucleic Acid Detection Kit	Nasopharyngeal swab/Blister fluids	50 tests/kit	– (20±5) °C	12 months	Available	RUO
74	R253T050B0C0	Influenza A Virus N4 Gene Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
75	R254T050B0C0	Influenza A Virus N5 Gene Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
76	R261T050B0C0	Influenza A Virus H12 Subtype Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
77	R263T050B0C0	Influenza A Virus H15 Subtype Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
78	R265T050B0C0	Influenza A Virus (H1N3) nucleic acid test kit		50 tests/kit	– (20±5) °C	9 months	N/A	RUO
79	R281T010B0C0	Rhinovirus Nucleic Acid Detection Kit	Sputum/Oropharyngeal swab	10/25/50 tests/kit	– (20±5) °C	-	N/A	RUO
80	R351T050B0C0	Influenza A H1 and Seasonal H3 Subtype Nucleic Acid Detection Kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
81	R26701T2048	Influenza A H1N1(2009) Nucleic Acid Test Kit	swabs and sputum	48Tests/Box	2-8 °C	24 months	Available	CE IVD

NO	CAT.NO	PRODUCT DESCRIPTION	Specimen	Kit Size	Storage Temperature	Shelf Life	English IFU	Certificate
82	R401T002B0C2	2019-Novel Coronavirus B.1.1.7(a) Mutant test kit	Sputum/ Sliva/ Nasopharyngeal swab/ Oropharyngeal swab /Bronchoalveolar lavage fluid/ Feces	50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
83	R401T003B0C2	2019-Novel Coronavirus B.1.351(β) Mutant test kit		50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
84	R401T004B0C2	2019-Novel Coronavirus P.1(y) Mutant test kit		50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
85	R401T008B0C2	2019-Novel Coronavirus B.1.1(y) 2(δ) Mutant test kit		50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
86	R401T005B0C2	2019-Novel Coronavirus B.1.617.2(δ) + B.1.617.1(k) Mutant test kit		50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
87	R401T006B0C2	2019-Novel Coronavirus .1.525(η) Mutant test kit		50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
88	R401T007B0C2	2019-Novel Coronavirus B.1.620 Mutant test kit		50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
89	R401T016B0C2	2019-Novel Coronavirus Omicron Mutant test kit		10/50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
90	R401T011B0C2	2019-Novel Coronavirus 8 Mutants test kit		10/50 tests/kit	– (20±5) °C	9 months	Available	CE IVD
<b>Tumor</b>								
91	T011T012B0C2	V600E mutation in B-raf gene test kit	FFPE	12tests/kit	– (20±5) °C	12 months	Available	NMPA/ CE IVD
92	T021T012B0C2	7 kinds of mutation in K-ras gene test kit	FFPE	12tests/kit	– (20±5) °C	12 months	Available	NMPA/ CE IVD
93	T071T012B0C0	EGFR 29 common mutation gene test kit	lung lesions, fresh frozen tissues, FFPE	12tests/kit	– (20±5) °C	12 months	Available	CE
<b>Sexual Transmitted Diseases</b>								
94	S015T024B0C2	<b>Human Papillomavirus nucleic acid amplification test kit</b> (Group 1: 6, 11, 42, 43, 44, 81 - FAM Group 2: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 68 - ROX Group 3: 26, 53, 66, 73, 82 - HEX/VIC Internal Control - Cy5)	Cervical Shedding Cells	24 tests/kit	– (20±5) °C Liquid	12 months	Available	NMPA/ CE IVD/ UKCA
	S015T048B0C2			48 tests/kit				
	S015T096B0C2			96 tests/kit				
95	S01502T2048	Internal Control - Cy5) - 1 rxn tube	Cervical Shedding Cells	48 tests/kit	<b>2-8°C</b> (Lyophilized)	24 months	Available	NMPA/ CE IVD/ UKCA
	S01502T2064			64 tests/kit				
	S01502T2096			96 tests/kit				
96	S014T024B0C2	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> (24 Genotypes: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 68, 26, 53, 66, 73, 82, 6, 11, 42, 43, 44 and 81) - 8 rxn tubes.	Cervical Shedding Cells	24 tests/kit	– (20±5) °C Liquid	12 months	Available	NMPA/ CE IVD/ UKCA
	S014T048B0C2			48 tests/kit				
97	S01402T2024		Cervical Shedding Cells	24 tests/kit	<b>2-8°C</b> (Lyophilized)	24 months	Available	
98	S021T024B0C2	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> (16, 18 and other high- and middle-risk Human Papillomavirus Types 16 - FAM 18 - HEX/VIC Other 16 pooled group - ROX Internal Control - Cy5)	Cervical Shedding Cells	24 tests/kit	– (20±5) °C Liquid	12 months	Available	NMPA/ CE IVD/ UKCA
	S021T048B0C2			48 tests/kit				
	S021T096B0C2			96 tests/kit				
99	S02104T2048	Internal Control - Cy5) - 1 rxn tube	Cervical Shedding Cells/Urine	24 tests/kit	<b>2-8°C</b> (Lyophilized)	24 months	Available	
	S02104T2064			48 tests/kit				
	S02104T2096			96 tests/kit				
100	S02701T2048	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> 16, 18/ 45, and other 11 pooled group (31, 33, 35, 39, 51, 52, 56, 58, 59, 66, 68) - 1 rxn tube	Cervical Shedding Cells/Urine	48 tests/kit	<b>2-8°C</b> (Lyophilized)	24 months	Available	CE IVD
	S02701T2096			96 tests/kit				
101	S02801T2048	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> (8 genotypes: 16, 18, 31, 33, 35, 45, 52 and 58) - 4 rxn tubes.	Cervical Shedding Cells	48 tests/kit	<b>2-8°C</b> (Lyophilized)	24 months	Available	CE IVD
102	S02901T2048	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> 16,18 and other 6 pooled group (31, 33, 35, 45, 52, 58) - 1 rxn tube	Cervical Shedding Cells	48 tests/kit	<b>2-8°C</b> (Lyophilized)	24 months	Available	CE IVD
103	S03001T2048	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> 16, 18, and other 16 pooled group (26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 73 and 82) - 1 rxn tube	Cervical Shedding Cells	48 tests/kit	<b>2-8°C</b> (Lyophilized)	12 months	Available	CE IVD
	S03001T2096			96 tests/kit				
104	S031T024B0C2	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> 16, 18, and other 12 pooled group (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68) - 1 rxn tube	Cervical Shedding Cells	24 tests/kit	– (20±5) °C Liquid	12 months	Available	CE IVD
	S031T048B0C2			48 tests/kit				
	S031T096B0C2			96 tests/kit				
105	S032T024B0C2	<b>Human Papillomavirus Nucleic Acid Amplification Test Kit</b> (31 Genotypes: 16, 18, 31, 43, 33, 35, 39, 73, 45, 51, 52, 56, 58, 59, 82, 68, 26, 53, 81, 6, 11, 42, 66, 40, 44, 54, 61, 62, 69, 70, 84) - 8 rxn tubes.	Cervical Shedding Cells	24 tests/kit	– (20±5) °C Liquid	12 months	Available	CE IVD
106	S073T050B0C0	<b>Herpes simplex virus types I and II typing kit</b>	Herpes fluid	50 tests/kit	– (20±5) °C	12 months	Available	CE IVD
107	S074T050B0C0	Herpes Simplex Virus Type II (HSV-II) (Quantitative)	Whole Blood	50 tests/kit	– (20±5) °C	12 months	Available	RUO
108	S075T050B0C0	Herpes Simplex Virus Types I Test Kit	whole blood	50 tests/kit	– (20±5) °C	12 months	Available	RUO
109	S081T050B0C0	Human parvovirus B19 nucleic acid test kit	Serum/Plasma	50 tests/kit	– (20±5) °C	12 months	Available	RUO
110	S08201T2024	<b>Urological and Reproductive Tract Microbial Nucleic Acid Test Kit</b> (31 pathogens, 8 reaction tubes)	Female cervical swabs, cervical exfoliated cell samples, and male urethral swabs	24 Tests/Kit	2-8 °C	12 months	Available	CE IVD
<b>Bacterial Infection &amp; Fungal Infection</b>								
111	B051T050B0C2	<i>Helicobacter pylori</i> Nucleic Acid Detection	gastric mucosal biopsy	50 tests/kit	– (20±5) °C	12 months	Available	RUO
112	B052T050B0C2	Clarithromycin, Quinolone-resistant and <i>Helicobacter pylori</i> nucleic acid amplification test kit	Gastric mucosal biopsies	50 tests/kit	– (20±5) °C	12 months	Available	NMPA/ CE IVD UKCA
113	B05201T2048	<b>Clarithromycin, Quinolone-resistant and <i>Helicobacter pylori</i> nucleic acid amplification test kit</b>	Feces	48 tests/kit	<b>2-8°C</b> (Lyophilized)	12 months	Available	CE IVD UKCA
114	B053T025B0C2	Human CYP2C19 Gene Polymorphism Detection Kit	WB/human tissue	25 tests/kit	– (20±5) °C	12 months	Available	RUO
115	B054T050B0C0	<i>Helicobacter pylori</i> Typing Nucleic Acid Detection Kit	gastric mucosal biopsy	50 tests/kit	– (20±5) °C	12 months	Available	RUO

NO	CAT.NO	PRODUCT DESCRIPTION	Specimen	Kit Size	Storage Temperature	Shelf Life	English IFU	Certificate
116	B055T010B0C2	<b>Multi-Antibiotic Resistant H.pylori Strain Real Time PCR Detection Kit (Furazolidone, tetracycline and amoxicillin-resistant)</b>	gastric mucosal biopsy	10 tests/kit	– (20±5) °C	12 months	Available	CE IVD
117	B081T050B0C0	<b>Group B Streptococcus (Streptococcus agalactiae) nucleic acid test kit</b>	Vaginal secretion/ Anal swab	50 tests/kit	– (20±5) °C	12 months	Available	CE IVD
118	B091T024B0C0	<i>Bordetella pertussis</i> nucleic acid test kit		24 tests/kit	– (20±5) °C	12 months	Available	RUO
119	B092T050B0C0	<i>Bordetella pertussis</i> and <i>Bordetella parapertussis</i> nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
120	B093T050B0C0	<i>Bordetella pertussis</i> , <i>B. parapertussis</i> , and <i>B. holmesii</i> nucleic acid test kit	Oropharyngeal swab	50 tests/kit	– (20±5) °C	12 months	Available	RUO
121	B161T050B0C0	Enterotoxigenic <i>Escherichia coli</i> nucleic acid test kit		50 tests/kit	– (20±5) °C	12 months	Available	RUO
122	B162T050B0C0	Enteroinvasive <i>Escherichia coli</i> nucleic acid test kit		50 tests/kit	– (20±5) °C	12 months	Available	RUO
123	B163T050B0C0	Enteropathogenic <i>Escherichia coli</i> nucleic acid test kit		50 tests/kit	– (20±5) °C	12 months	Available	RUO
124	B164T050B0C0	Enterohemorrhagic <i>Escherichia coli</i> test kit		50 tests/kit	– (20±5) °C	12 months	Available	RUO
125	B165T050B0C0	Enterogastric <i>Escherichia coli</i> typing kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	12 months	Available	RUO
126	B166T050B0C0	<b>Multiplex 5 Diarrheagenic Escherichia Coli Nucleic Acid test kit (Enterotoxigenic, Enteroinvasive, Enteropathogenic, Enterohemorrhagic and Enterogastric <i>Escherichia coli</i>)</b>		50 tests/kit	– (20±5) °C	12 months	Available	RUO
127	B173T050B0C0	<i>Vibrio Cholerae (group O139)</i>		50 tests/kit	– (20±5) °C	12 months	Available	RUO
128	B174T050B0C0	<i>Vibrio Cholerae (group O1)</i>		50 tests/kit	– (20±5) °C	12 months	Available	RUO
129	B175T050B0C0	<i>Vibrio Cholerae ctx Virulence Gene</i>		50 tests/kit	– (20±5) °C	12 months	Available	RUO
130	B176T050B0C0	<i>Vibrio Parahaemolyticus</i> and <i>Virulence Genes</i>	Feces/Anal swab	50 tests/kit	– (20±5) °C	12 months	Available	RUO
131	B177T050B0C0	<i>Vibrio Cholerae O1/O139</i> and <i>ctx Virulence Gene</i>		50 tests/kit	– (20±5) °C	12 months	Available	RUO
132	B181T050B0C0	<i>Salmonella</i> and <i>Shigella</i> nucleic acid test kit		50 tests/kit	– (20±5) °C	12 months	Available	RUO
133	B191T050B0C0	<i>Streptococcus pneumoniae</i> and Haemophilus influenza		50 tests/kit	– (20±5) °C	-	N/A	RUO
134	B201T050B0C0	<i>Legionella</i> nucleic acid test kit	Nasopharyngeal swabs/ Sputum/ Pleural effusion/ Lung tissue biopsy specimens	50 tests/kit	– (20±5) °C	12 months	Available	RUO
135	B221T050B0C0	<i>Haemophilus influenzae</i> type b nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
136	B231T050B0C0	<i>Scarlet streptococcus</i> nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
137	B241T010B0C0	<i>Neisseria meningitidis</i> Serogroup A		10 tests/kit	– (20±5) °C	12 months	Available	RUO
138	B242T010B0C0	<i>Neisseria meningitidis</i> Serogroup W135	whole blood /cerebrospinal fluid	10 tests/kit	– (20±5) °C	12 months	Available	RUO
139	B243T010B0C0	<i>Neisseria meningitidis</i>		10 tests/kit	– (20±5) °C	12 months	Available	RUO
140	B251T050B0C0	<i>Plesiomonas Shigelloides</i> nucleic acid test kit	feces or anal swab	50 tests/kit	– (20±5) °C	12 months	Available	RUO
141	B261T050B0C0	<i>Bacillus paratyphosus</i> A nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
142	B262T050B0C0	<i>Bacillus paratyphosus</i> B nucleic acid test kit		50 tests/kit	– (20±5) °C	-	N/A	RUO
143	B281T050B0C0	<i>Clostridium difficile</i> nucleic acid test kit	Feces/Anal swab	50 tests/kit	– (20±5) °C	12 months	Available	CE IVD
144	B28201T2048	<i>Clostridium Difficile</i> TcdA Gene And TcdB Gene Nucleic Acid Test Kit	Feces and anal swab	48 Tests/Box	2-8°C (Lyophilized)	24 months	Available	CE IVD/UKCA
145	B301T050B0C0	<i>Aeromonas hydrophila</i> nucleic acid test kit	feces and anal swab	50 tests/kit	– (20±5) °C	12 months	Available	RUO
146	B701T050B0C1	<i>Candida Auris</i> Nucleic Acid Test Kit	Dermal or environmental swabs	50 Tests/Box	– (20±5) °C	12 months	Available	CE IVD
147	B70201T2048	<i>Candida auris</i> and Antibiotic-Resistant Strain Nucleic Acid Test Kit	whole blood/ serum, plasma/swab	48 Tests/Box	2-8 °C	12 months	Available	CE IVD/UKCA
148	sample is ready	<b>Gastricintestinal infection Panel 1 (24 pathogens)</b> Viruses (6): rotavirus (groups A and B), norovirus (GI and GII), sapovirus (types I, II, IV and V), astrovirus, adenovirus group F (40 and 41), enterovirus Bacteria (13): enteroaggregative <i>Escherichia coli</i> (EAEC), enterotoxigenic <i>Escherichia coli</i> (ETEC), enteropathogenic/enterohemorrhagic <i>Escherichia coli</i> (EPEC/EHEC), <i>Shigella</i> /enteroinvasive <i>Escherichia coli</i> (EIEC), <i>Salmonella</i> , <i>Campylobacter</i> ( <i>Campylobacter jejuni</i> , <i>Campylobacter coli</i> / <i>Campylobacter uprasa</i> ), <i>Campylobacter jejuni</i> , <i>Yersinia enterocolitica</i> , <i>Vibrio cholerae</i> , <i>Aeromonas</i> ( <i>Aeromonas hydrophila</i> / <i>Aeromonas caviae</i> , <i>Aeromonas dakar</i> / <i>Aeromonas vermiculata</i> ), <i>Aeromonas hydrophila</i> , <i>Clostridium difficile</i> (toxin A and toxin B) <i>Staphylococcus aureus</i> Parasites (4): <i>Cryptosporidium</i> , <i>Enterosporidiosis</i> , <i>Entamoeba histolytica</i> , <i>Giardia lamblia</i> Fungi (1): <i>Candida albicans</i>	Human feces/anal swabs/vomit	24 tests/kit	– (20±5) °C	12 months	Available	CE IVD
<b>Others</b>								
149	M011T024A0C2	Human leukocyte antigen B27 nucleic acid test Kit	Whole Blood	24 tests/kit	– (20±5) °C	12 months	Available	NMPA/CE IVD
150	M09202T2048	Dengue virus nucleic acid test kit(with internal control)	Serum	48 Tests/Box	2-8°C (Lyophilized)	24 months	Available	CE IVD
151	M09601T2048	Dengue Virus Genotyping Nucleic Acid Test Kit	serum	48 Tests/Box	2-8°C (Lyophilized)	24 months	Available	CE IVD
152	M101T010B0C0	Zika Virus Nucleic Acid Test Kit	Serum/Plasma	10 tests/kit	– (20±5) °C	9 months	Available	RUO
153	M141T010B0C0	Rift Valley Fever Virus Nucleic Acid Detection Kit	Serum	10 tests/kit	– (20±5) °C	9 months	Available	RUO
154	M241T010B0C0	Chikungunya Virus Nucleic Acid Detection Kit	Serum	10 tests/kit	– (20±5) °C	9 months	Available	RUO
155	M10201T2048	Dengue/Zika/Chikungunya Triplex Nucleic acid test kit	Serum	48 Tests/Box	2-8°C (Lyophilized)	24 months	Available	CE IVD
156	M251T010B0C0	Yellow Fever Nucleic Acid Detection Kit	Serum	10 tests/kit	– (20±5) °C	9 months	Available	RUO
157	R60101T2050	Monkeypox Virus Nucleic Acid Detection Kit	Plasma/Rash Exudate	50 tests/kit	– (20±5) °C	12 months	Available	CE IVD/UKCA
158	TP060101T2024	Human Telomere Length Real Time PCR Detection Kit	Whole Blood/Throat Swab	24 tests/kit	– (20±5) °C	12 months	Available	CE IVD
							Date:	27-Jun-25

**JIANGSU MOLE BIOSCIENCE CO., LTD**

Tel: 0571-87209310 Fax:0523-82205058

info@molechina.com

www.molechina.com

